

Product Manual

Ultra-Helium Dewar[™]Cryogenic Tanks 60, 100, 250 & 500 Liter Models



Designed and Built by:

Chart Inc.

1300 Airport Dr. Ball Ground, GA 30107 USA (800) 400-4683



Contents

| Revision Log | iv |
|---|----|
| Preface | .1 |
| General | |
| Product Highlights. | .1 |
| Product Manual | |
| Terms | |
| Acronyms / Abbreviations. | |
| | 2 |
| Safety | |
| | |
| General | |
| Safety Summary. | |
| Oxygen Deficient Atmospheres | |
| Oxygen Enriched Atmospheres | |
| | |
| Handling | |
| Handling | .7 |
| Transportation by Truck | .8 |
| Transporting by Air Cargo | .8 |
| Equipment Description | 0 |
| General | |
| Thermal Insulation | |
| Plumbing Components | |
| Optional Equipment | |
| Helium Transfer Line. | |
| Standard Features. | |
| Special Features | |
| | |
| Operation | |
| General | |
| Liquid Nitrogen Pre-Cool and Helium Purge | |
| Initial Liquid Helium Filling | |
| Liquid Helium Filling | |
| Liquid Withdrawal. | |
| External Pressurization | |
| Product Loss from Liquid Transfer | |
| Storage. | |
| Liquid Level Reading | |
| Troubleshooting | 18 |
| Specifications | 19 |
| Drawings and Parts Lists | |
| Ultra-Helium Dewar 60 / 100 | 20 |
| 60 / 100 Parts List | 21 |
| Ultra-Helium Dewar 250 / 500 | |
| 250 / 500 Parts List | 23 |
| Wiring Diagram. | 24 |
| Pressure Builder | 25 |

Revision Log

| Revision Level | Date | Description |
|-----------------------|------------|-----------------------|
| A | 10/01/2004 | Original |
| В | 12/15/2016 | Updated to new format |
| | | |
| | | |
| | | |



Preface

General

The Ultra-Helium Dewar cryogenic tanks are designed and built for reliable transport. They are light, maneuverable and durable, while providing superior thermal performance. The unique neck tube design provides proven support during transportation. The outboard caster base provides maximum stability in a compact design.

Available in sizes ranging from 60 to 500 liters. The Ultra-Helium Dewar tanks are suitable for air transport (IATA conforming) with the optional absolute pressure relief valve. All models are 100% non-magnetic for Magnetic Resonance Imagery (MRI) service.

Product Highlights

The Ultra-Helium Dewar cryogenic tank controls are conveniently located on the top, with nesting fill couplings to accept various standard transfer line sizes. The optional electric pressure builder can quickly increase pressure for liquid transfer while maintaining low heat leak. It has two preset ranges (4 or 8 psig/0.3 or 0.6 barg) for efficient liquid helium withdrawal.

- Maximum durability and lightweight
- Outstanding thermal performance
- Large ball valves for up to 3/4" (19mm) transfer lines

Product Manual

The Ultra-Helium Dewar Product Manual is designed to be used in conjunction with Ultra Dewar 60, 100, 250, and 500 liter liquid helium vessels. It should be thoroughly read and understood by anyone that operates, or is exposed to this equipment. If there are any questions regarding the operation of the cylinder, contact Chart's Technical Service division at 1-800-400-4683.

The safety requirements for handling and transporting liquid helium are shown in the Safety section. It is imperative that all persons having contact with the Ultra-Helium Dewar become thoroughly familiar with all maintenance, safety precautions, and procedures contained in this product manual.

The Handling section will provide the operator with the proper handling and transporting instructions.

The Equipment Description section will familiarize the operator with the working components of the various models in the Ultra-Helium Dewar series.

The Operation section contains detailed information on the general operation of the Ultra-Helium Dewar container including filling, liquid withdrawal and storage. A Troubleshooting section is also included in this section.

Please refer to the Specifications section for a complete listing of part numbers, drawings and other technical information.

Terms

Throughout this manual safety precautions will be designated as follows:



Warning! Description of a condition that can result in personal injury or death.



Caution! Description of a condition that can result in equipment or component damage.



: A statement that contains information that is important enough to emphasize or repeat.

Acronyms / Abbreviations

The following acronyms / abbreviations are used throughout this manual:

| BAR | Pressure (Metric) |
|-------|---|
| BARG | Pressure (Metric) Gauge |
| BTU | British Thermal Unit |
| DOT | Department of Transportation |
| G | Grams |
| He | Helium |
| IATA | International Air Transport Association |
| Kg | Kilogram |
| Kj/Kg | Kilojoule/Kilogram |
| MAWP | Maximum Allowable Working Pressure |
| MRI | Magnetic Resonance Imagery |
| PB | Pressure Builder |
| PPM | Parts Per Million |
| PSI | Pounds per Square Inch |
| PSIG | Pounds per Square Inch (Gauge) |
| | |



Safety

Introduction to Helium

Liquid Helium has unique properties that require the usage of this specially designed Ultra-Helium Dewar cryogenic container. This section of the manual will discuss the properties of helium and how the Ultra-Helium Dewar cryogenic equipment is affected by them

Helium is found in the atmosphere in only small trace amounts (5 ppm) and is not economically obtained by air separation. Helium also exists in larger amounts (2%) in certain types of natural gas wells.

Helium is chemically inert. It has no color, odor or taste. Helium is non-flammable and only slightly soluble in water. Helium is present in two stable isotopes; Helium 3 (He³) and Helium 4 (He⁴). He⁴ is the most common, and makes up over 99% of the helium being used today. This manual will discuss only the properties of He⁴ and will generically refer to it as helium.

Liquid helium exists at very low temperatures (-453°F, -269°C) at atmospheric pressures, a temperature that is only a few degrees above absolute zero. Helium will not solidify as the temperature approaches absolute zero, therefore making helium an ideal liquid for low temperature physics. The extremely low temperature requires that additional safety precautions be observed while operating the equipment. Liquid helium should not be allowed to come in contact with air. The extreme cold will liquefy the oxygen in the air and cause a safety hazard associated with a high oxygen environment.

Liquid helium has a very low latent heat of vaporization (8.8 BTU/lb, 20.5 kj/kg). This is 1/110 that of water. Liquid helium boils away very easily and requires special equipment and procedures to reduce the heat input that causes evaporation. The Ultra-Helium Dewar container takes advantage of helium's high sensible heat (640 BTU/lb, 1,448.6 kj/kg from boiling point to 70°F, 21°C) to offset its low latent heat. The Ultra-Helium Dewar container insulation system has radiation shields that intercepts the incoming heat and transfers it into the cold gas that is venting out of the container.

Liquid helium has a very low critical pressure (18 psig/1.2 barg) where it exists as both liquid and gas and requires very low operating pressures. The container must be a closed system and have a pressure greater than atmospheric. If the containers valves were left open and the atmospheric pressure became greater, air would rush into the vessel and freeze. The frozen water vapor in the air would quickly plug the vessel openings and create a dangerous situation. The

Ultra-Helium Dewar protects itself with a series of relief valves that are attached to various parts of the plumbing. It is extremely important to make sure all valves are closed and that all relief valves close with a tight seal.

General

This section of the manual deals with the safety precautions that are necessary with cryogenic equipment. The potential hazards in handling liquid helium stems mainly from the physical properties:

- The liquid is extremely cold (helium is the coldest of all cryogenic liquids)
- The ultra-low temperatures of liquid helium will condense and solidify air
- Very small amounts of liquid helium are converted into large amounts of helium gas
- · Helium gas is non-life supporting

While Chart equipment is designed and built to rigid standards, no piece of mechanical equipment can ever be made 100% safe. Strict compliance with proper safety and handling practices are necessary when using a liquid cylinder or other compressed gas equipment. We recommend that all of our customers reemphasize safety and safe handling practices to all their employees and customers. While safety features have been designed into the unit and safe operations are anticipated, it is essential that the user of these liquid cylinders carefully read to fully understand all WARNINGS and CAUTIONS listed in this safety section and enumerated below.



Warning! Keep air and other gases away from liquid helium. The low temperatures of liquid helium or cold gaseous helium can solidify any other gas. Solidified gases or liquid that is allowed to form on the plumbing can plug pressure relief passages and relief valves. Plugged passages are hazardous because of the continual need to relieve excess pressure. Always store and handle liquid helium under positive pressure. Use liquid helium in a closed system to prevent the infiltration of air or other gases.



Warning! Keep exterior surfaces clean to prevent combustion. Atmospheric air will condense on exposed helium piping. The nitrogen, having a lower boiling point than the oxygen, will evaporate first from the condensed air, leaving an oxygen enriched liquid. This liquid may drip or flow to nearby surfaces. These surfaces must be clean to "oxygen clean" standards to prevent a possible combustion.

Safety Summary

In the event a liquid cylinder is inadvertently dropped, tipped over, or abused, slowly raise it to its normal vertical position. Immediately open the vent valve to release any excess pressure in a safe manner. As soon as possible, remove the liquid product from the vessel in a safe manner. If the vessel has been used in oxygen service, purge it with an inert gas (nitrogen). If damage is evident or suspected, return to Chart prominently marked "LIQUID CYLINDER DROPPED, INSPECT FOR DAMAGE."



Caution! Before removing cylinder parts or loosening fittings, completely empty the liquid cylinder of liquid and release the entire vapor pressure in a safe manner. External valves and fittings can become extremely cold and may cause painful burns to personnel unless properly protected.

Oxygen Deficient Atmospheres



Warning! Keep the equipment area well ventilated. Although helium is non-toxic it can cause asphyxiation in a confined area without ventilation. Any atmosphere that does not contain enough oxygen for breathing can cause dizziness. unconsciousness, or even death. Helium, being colorless, odorless, and tasteless, cannot be detected by human senses. Without adequate ventilation, helium will displace the air and give no warning that a nonlife supporting atmosphere is present. Store liquid helium in a well ventilated area.

The normal oxygen content of air is approximately 21%. Depletion of the oxygen content in air, either by combustion or by displacement with inert gas, is a potential hazard and users should exercise suitable precautions.

One aspect of this possible hazard is the response of humans when exposed to an atmosphere containing only 8 to 12% oxygen. In this environment unconsciousness can be immediate with virtually no warning.

When the oxygen content of air is reduced to about 15% to 16%, the flame of ordinary combustible materials, including those commonly used as fuel for heat or light, may be extinguished. Somewhat below this concentration, an individual breathing the air is mentally incapable of diagnosing the situation because the onset of symptoms such as sleepiness, fatigue, lassitude, loss of coordination, errors in judgment and confusion can be masked by a state of "euphoria," leaving the victim with a false sense of security and well being.

Human exposure to atmosphere containing 12% or less oxygen leads to rapid unconsciousness. Unconsciousness can occur so rapidly that the user is rendered essentially helpless. This can occur if the condition is reached by an immediate change of environment, or through the gradual depletion of oxygen.

Most individuals working in or around oxygen deficient atmospheres rely on the "buddy system" for protection obviously the "buddy" is equally susceptible to asphyxiation if he or she enters the area to assist the unconscious partner unless equipped with a portable air supply. Best protection is obtained by equipping all individuals with a portable supply of respirable air. Life lines are acceptable only if the area is essentially free of obstructions and individuals can assist one another without constraint.

If an oxygen deficient atmosphere is suspected or known to exist:

- 1. Use the "buddy system." Use more than one "buddy" if necessary to move a fellow worker in an emergency.
- 2. Both the worker and "buddy" should be equipped with self-contained or airline breathing equipment.

Oxygen Enriched Atmospheres



Warning! In oxygen enriched atmospheres flammable items burn vigorously and could explode. Certain items considered non-combustible in air may burn rapidly in such an environment.

An oxygen enriched atmosphere occurs whenever the normal oxygen content of air is allowed to rise above 23%. While oxygen is nonflammable, ignition of combustible materials can occur more readily in an oxygen rich atmosphere than in air; and combustion proceeds at a faster rate although no more heat is released.

It is important to locate an oxygen system in a well ventilated location since oxygen rich atmospheres may collect temporarily in confined areas during the functioning of a safety relief device or leakage from the system.

Oxygen system components, including but not limited to, containers, valves, valve seats, lubricants, fittings, gaskets and interconnecting equipment including hoses, shall have adequate compatibility with oxygen under the conditions of temperature and pressure to which the components may be exposed in the containment and use of oxygen. Easily ignitable materials shall be avoided unless they are parts of equipment or systems that are approved, listed, or proven suitable by tests or by past experience.

Compatibility involves both combustibility and ease of ignition. Materials that burn in air may burn violently in pure oxygen at normal pressure, and explosively in pressurized oxygen. In addition, many materials that do not burn in air may do so in pure oxygen, particularly when under pressure. Metals for containers and piping must be carefully selected, depending on service conditions. The various steels are acceptable for many applications, but some service conditions may call for other materials (usually copper or its alloy) because of their greater resistance to ignition and lower rate of combustion.

Similarly, materials that can be ignited in air have lower ignition energies in oxygen. Many such materials may be ignited by friction at a valve seat or stem packing, or by adiabatic compression produced when oxygen at high pressure is rapidly introduced into a system initially at low pressure.

Personal Protective Equipment (PPE)



Warning! Cover eyes and exposed skin. Accidental contact of liquid helium or cold issuing gas with the eyes or skin may cause freezing injury similar to a burn. Protect the eyes and cover the skin where possible contact with cold liquids or gas exists.

The following personal protective equipment is recommended when working around cryogenic liquid:

- Safety glasses with side shields to prevent cryogenic liquid from splashing into the eyes
- Chemical / Liquid resistant gloves to prevent cryogenic burns on exposed hands
- Long sleeve shirts to protect the arms
- Cuffless trousers worn over closed shoes



Handling

Handling

The Ultra-Helium Dewar cryogenic tanks are designed to be portable and are equipped with square roller bases and handles. The various models have three different styles of portable handling systems.



Warning! DO NOT handle these containers in any other fashion than what is described in this manual. DO NOT transport or store the containers on their side. DO NOT attempt to lift the container by any other means than described in this manual. Failure to observe proper handling procedures could cause the container to tip over causing possible injury.



Caution! The Ultra-Helium Dewar containers are designed to be moved over clean, flat, hard surfaces only. DO NOT attempt to roll the container over aberrations in the floors surface that could lock a wheel and tip the container.

The Ultra-Helium Dewar 60 and 100 models have square bases with two fixed, and two swivel casters. These models have a handle mounted on each side of the outer container. Figure 1 illustrates the handles, base, and casters for the Ultra-Helium Dewar 60 and 100 models.



The Ultra-Helium Dewar 250 model also has a square base with two fixed, and two swivel casters. These models are equipped with a pull handle assembly (see Figure 2).





The Ultra-Helium Dewar 500 model has four bracketmounted casters equally spaced around the circumference of the vessel. It also has two fixed, two swivel casters and four side mounted grab handles (see Figure 3).



Figure 3

Figure 1

Transportation by Truck

The Ultra-Helium Dewar cryogenic containers are low pressure (below 25 psi) portable containers that are not required to meet DOT requirements for the transport of liquid helium. The container does require the proper labeling for low temperature liquid helium and should be handled carefully and with the respect due any cryogenic container. When transporting by truck:

- Use a lift gate to load and unload the container from the truck
- Stand next to the container and steady it while it is on the lift gate
- Use nylon straps to secure the container in the truck



Caution! DO NOT use chains to secure the containers. The chain will damage the outer vacuum jacket.

Transporting by Air Cargo

The Ultra-Helium Dewar cryogenic containers are designed to be transported by air cargo and meet the requirements of International Air Transport Association (IATA). When an air shipment is required the shipper should make advanced arrangements with the cargo manager. This should avoid any questions or unnecessary delays when the container is presented for shipment.

The Ultra-Helium Dewar container must be secured, or packed in a proper shipping crate or frame, to prevent it from tipping or moving during transport. The container should have the following shipping notices on it, at 120° intervals, in clear view of the handlers:

DO NOT DROP

THIS SIDE UP (with arrows)

CARGO AIRCRAFT ONLY

HELIUM, LIQUID, LOW PRESSURE (and the proper DOT product labeling)

The change of atmospheric pressure, even in a pressurized compartment, can cause a problem with the relief valves. Ordinary relief valves can freeze and plug up following rapid ejection of cold gas following altitude changes. An absolute pressure relief valve should be used for all air transportation of liquid helium. The relief valves should be examined upon arrival for normal operation to assure that ice plugs have not occurred.

CHART

Equipment Description

General

The Ultra-Helium Dewar cryogenic containers are designed for the specific application of transporting and storage of liquid helium. These containers should not be used for other cryogenic liquids. They do not require liquid nitrogen shielding to achieve good thermal performance. They are constructed with a stainless steel inner and outer vessel. The containers are attached to roller bases that make them easily transported. Refer to the components illustration (Figure 4) for the nomenclature of the Ultra Dewar parts.





Thermal Insulation

The specially designed insulation system is composed of multiple layers of paper and foil that use a low vacuum to provide an extremely low heat leak between the inner and outer vessel. The system also employs radiation shields that take advantage of the sensible heat of the boiling liquid helium to reduce the normal boil-off of the product. As the helium boils it forces gas up the passage between the stainless steel neck tube to vent at the shield flow valve. The cold helium gas is raised in temperature by absorbing heat from the radiation shields as it travels up the neck tubes. The insulation and vacuum systems are designed to be selfmaintaining for years. The vacuum is sealed with an o-ring sealed plug that also acts as a relief device in case there is an inner vessel leak.



Caution! DO NOT tamper with the vacuum pump out plug. The rapid loss of container vacuum can cause damage to the insulation system.

Plumbing Components

The operational controls on the Ultra-Helium Dewar container are very simple. They are located on the top of the container and protected by the containers protection ring. The controls operate all of the functions of the Ultra-Helium Dewar containers operations. Refer to the schematic (Figure 5 on the next page) for the plumbing component locations.

The fill and withdrawal ball valve is located on the top of the tank in the vertical position with quick couplings attached to it. A transfer line can be passed through the coupling and ball valve and into the inner vessel of the container. The transfer line is then sealed in place with the quick coupling. Each container has a variety of quick couplings to allow for different sizes of transfer lines.

The vent valve is located so that it can be used to depressurize the tank during filling or pressurize the tank from another source for liquid transfer.



Caution! DO NOT leave the vent valve open after filling or transfer operation.

The relief valves are located so that if one stops operating from an ice plug the other will function properly. The primary relief valve is set at 10 psi. The alternate relief valve is set at 12 psi. Each one protects the inner pressure vessel, but from different paths.



Warning! DO NOT use relief valves that are set at higher pressure than allowed by the Ultra-Helium Dewar equipment. Refer to the Specifications section of this manual for the maximum allowable pressure setting.



Figure 5

The shield flow valve is set at 1 psi. It is located behind an isolation valve that allows it to be turned off while the container is being filled or pressurized for liquid transfer. This valve must be open for the insulation shields to be active.



The thermal performance of the Ultra-Helium Dewar will be reduced if the shield flow valve is closed or malfunctioning.

The container's pressure gauge is designed so that it will show vacuum as well as pressure. There are times when changes in the outside atmospheric pressure can cause the inner vessel to have a lower pressure than the atmospheric pressure.



Warning! DO NOT open any valves if the tank pressure reads below 0 on the pressure gauge. This would allow air to rush into the vessel and form an ice plug.

Optional Equipment

Helium Transfer Line

When using cryogenic fluids with a low heat of vaporization such as helium and hydrogen, it is necessary to have an efficient transfer apparatus. To meet these operating requirements, transfer tubes are carefully assembled using all welded stainless steel construction, and are super-insulated to provide the optimum in thermal performance.

All Chart transfer tubes are custom built to meet the customer's specific requirements for configuration, diameter, and length. Special features may be added as necessary, and fast delivery is available from our complete stock of components.

Standard Features

The inner configuration of all assemblies, whether they are rigid or flexible, is stainless steel of minimum wall thickness to reduce cool down losses and heat leak. Outer jackets are likewise stainless steel, capable of being evacuated to less than 10⁵ mm absolute pressure for minimum heat leak. Exceptions are taken when the customer requirements deviate from standard practices.

The evacuated space between the inner line and the outer jacket is super-insulated using aluminized mylar as a radiation barrier.

All tubes are helium mass spectrometer leak tested for a leak rate of less than 10⁶ micron cubic feet per hour, and are fully evacuated before shipping to ensure the unit is ready for immediate use upon receipt. An evacuation value is provided and can be used for occasional re-evacuation of the insulation space by the user.

Special Features (at extra cost)

- 1. **Bayonet Connectors.** Bayonet connectors can be provided for disassembly of the tube at any convenient location. These connectors can be easily assembled and disassembled in the laboratory, and do not require repumping of the insulating vacuum after each usage.
- 2. Flexible Sections. Chart can provide flexible sections on the horizontal legs of these tubes. These sections can be from 12" to 12' long as required by the customer. Bayonet couplings are also available with the flexible line sections.
- **3.** Vacuum Insulated Valves. Chart can provide 1/4" angle pattern super-insulated valves on the supply leg of these transfer tubes. This allows the on/off transfer type operation without loss of head pressure in the supply container.
- 4. Exchangeable Tips. Interchangeable uninsulated leg extensions with threaded connectors can be provided for versatile operation with different size containers and dewars.

When ordering, specify options required and specify dimensions A through E. Recommended combination of B/D and C/E to fit Ultra-Helium Dewar containers are as follows:

| | В | D | С | E |
|-----------------|-----|-----|-----|-----|
| Ultra-Dewar 60 | 29" | 13" | 29" | 13" |
| Ultra-Dewar 100 | 29" | 22" | 29" | 22" |
| Ultra-Dewar 250 | 33" | 30" | 33" | 30" |
| Ultra-Dewar 500 | 30" | 35" | 30" | 35" |



Figure 6

5. Optional Side Transfer Tube.



Figure 7



Operation

General

The Ultra-Helium Dewar cryogenic containers should be operated and handled as a piece of specialized equipment. The operational personnel should be thoroughly trained in the vessels operation and the nature of liquid helium. Rough or careless handling must be avoided. Review the Safety section before each operation to reinforce the proper cautions and warnings.



Warning! Liquid helium containers and lines must be cleaned and purged before use.

It is necessary to purge the inner vessel and its related plumbing of the nitrogen gas that is normally in a new tank. The transfer line must also be purged. Liquid helium would solidify any air in it and cause an ice plug.

The inner vessel must also be cooled down to near liquid helium temperatures before it will accept liquid helium. Filling with liquid helium will cool the vessel and it's insulation so that it will accept liquid helium. The vent gas can be recovered for re-liquefaction during this process. Cooling the vessel with liquid helium without vent gas recovery can be very costly.

Liquid nitrogen can be used in conjunction with the purging operation to cool the vessel to -320°F (-140°C) before liquid helium is introduced into the vessel. New dewars are at ambient temperature and should be pre-cooled to prevent large filling losses. Use the following procedure to purge and cool with liquid nitrogen.



Liquid Nitrogen Pre-Cool and Helium Purge

(Refer to Figure 8)

- 1. Insert a nitrogen fill stinger into the inner vessel until it reaches the bottom. The fill stinger should pass through the proper size quick coupling and the open fill valve (Item B).
- 2. Open the vent valve (Item A).
- 3. Close the shield flow isolation valve (Item C).
- 4. Fill the vessel to 10-20% of its full capacity by weight with liquid nitrogen.



: Vigorous venting will occur since the vessel is at ambient temperature.

- 5. Allow the liquid nitrogen to completely cool the inner vessel for at least 24 hours.
- 6. Attach a gaseous nitrogen pressure line to the vent valve (Item A) and pressure the vessel to 5 psi (0.34 bar).
- 7. Remove the liquid nitrogen through the fill stinger.
- 8. Examine the container for vacuum integrity. Check that the pump-out plug is in place to confirm vacuum.



Caution! Evacuation of the inner vessel of a container that has lost its vacuum insulation will cause the inner vessel to collapse and require a major repair.

- 9. Attach a vacuum pump to the vent valve (Item A).
- 10. Open the vent valve and close all other valves.
- 11. Evacuate the container and close the vent valve. Remove the pump.
- 12. Connect a dry helium gas source to the vent valve and break the vacuum back with helium gas.
- 13. Slowly raise the pressure to 1 psi (0.069 bar) to purge shield flow valve. Close shield flow isolation valve. Pressurize to 12 psi (0.82 bar) to purge the relief valves and vent system.
- 14. The vessel is now ready for filling with liquid helium.

Figure 8

Initial Liquid Helium Filling

The initial filling of the Ultra-Helium Dewar container is done in the same manner as described in the next paragraph on Liquid Helium Filling. It is very important to observe the following for the initial fill:

- 1. Make sure that the vessel was properly purged of air.
- 2. Make sure there is no residual liquid nitrogen if that was used for cool down.
- 3. Examine the vessel after filling.

Liquid Helium Filling

It is important that the contents of the Ultra-Helium Dewar container be measured as the filling operation is taking place. Liquid helium's heat of evaporation is so low that it is possible to go through the filling operation without collecting any liquid helium into the vessel. The heat input of any of the filling components or vessel could be high enough to vaporize the liquid before it reaches the vessel. The filling must be done through vacuum jacketed piping, valves, and a vacuum jacketed transfer tube. The best way to confirm the filling operation is to fill the Ultra-Dewar on a scale. The weight will start to increase as the vessel begins to accept liquid.

The filling should be done at low pressure. The optimum filling pressure varies with the configuration of each filling system. A balance between filling speeds and product loss will have to be maintained. The higher the pressure the faster the fill, but the more helium will be vaporized as the liquid helium depressurizes (see Figures 10 and 11). The average system filling pressure is between 2 and 8 psi (0.136 bar and 0.54 bar).

The following filling operation is done on a filling scale that is accurate to .1 pounds (45.4g) over the weight range of 100 to 1000 pounds (45.4 kg to 453.5 kg). This accuracy is necessary to indicate liquid helium, which weighs 0.28 lbs/ liter (.12 kg/liter). The operation is as follows (refer to Figure 8):

- 1. Purge and pre-cool the Ultra-Helium Dewar container as previously described.
- 2. Place the container on the filling scale. Record the empty weight.

- 3. Compare the empty weight with the tare weight that is stamped on the data plate. The difference could be cold helium gas (Cold helium gas at 0 psi equals 0.004 lbs/ liter or 0.17 kg/liter). If the difference is too large there could be residual helium in the vessel. If the vessel was just purged with liquid nitrogen, there could be residual nitrogen in the vessel. Re-purge the vessel if that is the case.
- 4. Purge the fill transfer tube with gaseous helium.
- 5. Open the fill valve (Item B). Insert the filling transfer tube through the appropriate quick connect (Item D). The filling tube should be pushed to the bottom of the inner vessel and tightened with the quick connect. Record the new weight. Subtract the vessel empty weight from the weight with the transfer hose. This will be the weight of the transfer hose. Add the transfer hose weight to the vessel tare weight and the desired amount of helium. This will be the target full weight.
- 6. Connect the vapor recovery system to the vent valve (Item A). If a vapor recovery system is not being used, the vent valve will be used to relieve helium-filling pressure. Venting will take place during the fill.
- 7. Open the vent valve (Item A) and the liquid helium transfer valve on the liquid source.
- 8. Watch the weight on the filling scale. Terminate the fill when the scale reads the target fill weight.
- 9. Close the transfer line valve. Allow the tank to vent down. Remove the filling tube and close the filling valve (Item B). Replace the plug in the quick connect that was used for filling. Close the vent valve (Item A).
- 10. Open the shield flow isolation valve (Item C). Make sure that all other valves are closed and that the tank is venting from the shield flow valve (Item A).

Liquid Withdrawal

Liquid helium can be withdrawn from the Ultra Dewar containers through the fill/withdrawal ball valve located on the top flange of the tank. The entire transfer must be made through vacuum insulated lines and into a cold vacuum insulated container. The most common procedure is to use a vacuum jacketed transfer line that can be provided from Chart (Technical Service 1-800-400-4863). Use the following procedure to transfer liquid helium from one tank to another without external pressurization (refer to Figures 8 and 9).



Figure 9

- 1. Close all valves on both the delivery and receiving tank.
- 2. Remove the quick coupling plugs on both tanks. Make sure that the proper size quick coupling (Item D) is being used for the transfer line.
- 3. Place the ends of the transfer line into the quick coupling on both tanks and insert them until they touch the fill/withdrawal ball valve (Item B). Tighten the quick couplings.



The boil-off, caused by the warm transfer line being inserted into the helium tank will cause immediate pressure rise and transfer of liquid helium. The transfer line should have a flow valve in it or be installed simultaneously into both tanks.

- 4. Open the ball valves (Item B) on both tanks and insert the transfer line into the tanks. Tighten the quick couplings (Item D).
- 5. Open the vent valve (Item A) on the receiving tank. Liquid transfer will take place as helium gas is vented from the receiving tank.
- 6. The liquid transfer can be measured by using a scale as described in the Liquid Helium Filling section of this chapter.

When the transfer is complete:

- Remove the transfer line until it just clears the fill/ withdrawal valve. Close the fill/withdrawal valve (Item B) and remove the line completely.
- 8. Close the receiving tanks vent valve. Open the shield flow valve (Item C) and let the tank vent to 1 psi (0.07 bar).
- 9. Install the quick coupling plugs and make sure that all valves except the shield flow valve are closed.

External Pressurization

Liquid transfer rate is a function of pressure. The higher the pressure the faster the transfer will occur. However, liquid helium gives off large quantities of gas as it is depressurized. The higher the transfer pressure the more liquid will be lost as gas. A balance is needed between transfer speed and product loss.

Pressure can be raised on the Ultra-Helium Dewar containers by connecting a high-pressure helium gas bottle with regulator to the Ultra-Helium Dewar's vent connection. The pressure should be adjusted between 4 and 8 psi (0.27 bar and 0.54 bar). The tank's pressure gauge and vent valve should be used to control the pressure. Pressures that exceed 10 psi will activate the tanks relief valves and cause large product losses.

An electronic pressure building (PB) system is available for some models of Ultra-Helium Dewar containers. The electric PB is designed to raise and maintain the pressure in the container to either 4 or 8 psi (0.27 bar or 0.54 bar). It indicates when the power and the pressure-building heater are on. The pressure builder should not be turned on until the transfer line is completely installed in both tanks (Step 4 in the withdrawal procedure).



Caution! The pressure builder must be turned off when the transfer is complete. The heater can be burned out if it is left on after the transfer.



Figure 10

Product Loss from Liquid Transfer

The Ultra-Helium Dewar container is a specially insulated tank that is designed to have very low product loss due to the normal heat leak into the vessel. There are other product loss factors that must be considered when making a liquid transfer. The introduction of the warm transfer line into the cold liquid helium will vaporize a certain amount of product. The cool down of the transfer line and receiving tank will use a certain amount of product. The addition of pressure to the tank from an outside source or the pressure builder will cause a certain amount of liquid to boil and increase the saturation pressure of the remaining liquid.

Flash loss due to pressure drop through the transfer line may be estimated by use of Figure 10 "Helium Flash Rate". Depressurization loss of liquid in the container may be estimated by use of Figure 11 "Helium Depressurization Losses".

When the supply of gaseous product is the primary operation of the tank, external vaporizers and an additional regulator must be added after the gas use valve (HCV-13) to step down the pressure to the gas application. This regulator (PCV-5) is found in the final line option provided by Chart.

For example: Assume a helium container is discharging at a constant pressure of 5 psig (0.34 bar). From Figure 10 the flash loss is approximately 13.8% of the liquid entering the transfer tube. From Figure 11 the loss from depressurizing the container is approximately 12.5% of the liquid remaining in the container.

For best transfer efficiency, the withdrawal should be started and maintained with as low a pressure as practical. Too low a pressure will require a longer time to make a transfer and thus permit heat leak in the transfer system to become excessive. A balance between effects of heat leak and depressurization may generally be attained by operating in a pressure range of 2 to 8 psig (0.13 to 0.54 barg).



Figure 11

Storage

Maximum storage efficiency is obtained when the liquid helium is maintained near atmospheric pressure. If liquid helium is stored in a closed vessel, the normal heat leak will cause a pressure rise, which ultimately must be relieved with accompanying blow down losses (refer to Figure 10)

This container is designed so that the refrigeration available in the boil-off gas from the liquid helium is used to cool multiple heat sink temperature barriers, which intercept the heat in-leak toward the helium reservoir. For operation, keep the 1 psi relief valve open and close all other vent valves. This forces the boil-off gas to circulate through the space between the necks. The relief valve is set to prevent air and/ or moisture from entering the vessel while permitting the boil-off gas to escape.

Liquid Level Reading

The most accurate means of determining the amount of helium in the Ultra Dewar is to weigh the container. The difference between the recorded weight and the tare weight is the amount of helium. The level of helium in the container can also be measured with a dipstick that is inserted through the fill ball valve to the surface of the liquid. The difference between the measured dipstick and the total length of the vessel can be compared to the contents chart for the associated volume of helium.

| Calibration Chart | | | | | | |
|--------------------|----------------------------|-----------------------------|-----------------------------|-----------------------------|--|--|
| Inches from Bottom | Ultra-Dewar 60 (liters) | Ultra-Dewar 100 (liters) | Ultra-Dewar 250 (liters) | Ultra-Dewar 500 (liters) | | |
| 0 | 0 | 0 | 0 | 0 | | |
| 1 | 1 | 1 | 2 | 2 | | |
| 2 | 4 | 4 | 7 | 7 | | |
| 3 | 9 | 9 | 15 | 15 | | |
| 4 | 14 | 14 | 24 | 26 | | |
| 5 | 19 | 19 | 34 | 39 | | |
| 6 | 24 | 24 | 44 | 54 | | |
| 7 | 29 | 29 | 55 | 70 | | |
| 8 | 34 | 35 | 66 | 87 | | |
| 9 | 39 | 40 | 76 | 104 | | |
| 10 | 44 | 45 | 87 | 122 | | |
| 11 | 50 | 50 | 98 | 140 | | |
| 12 | 55 | 55 | 109 | 157 | | |
| 13 | 60 | 60 | 119 | 175 | | |
| 14 | 65 | 65 | 130 | 192 | | |
| 15 | 68 | 71 | 141 | 210 | | |
| 16 | 69 | 76 | 152 | 228 | | |
| 17 | | 81 | 162 | 245 | | |
| 18 | | 86 | 173 | 263 | | |
| 19 | | 91 | 184 | 280 | | |
| 20 | | 96 | 195 | 298 | | |
| 21 | | 101 | 205 | 316 | | |
| 22 | | 106 | 216 | 333 | | |
| 23 | | 109 | 227 | 351 | | |
| 24 | | 110 | 242 | 369 | | |
| 25 | | | 248 | 386 | | |
| 26 | | | 258 | 404 | | |
| 27 | | | 266 | 421 | | |
| 28 | | | 272 | 439 | | |
| 29 | | | 275 | 457 | | |
| 30 | | | | 474 | | |
| 31 | | | | 490 | | |
| 32 | | | | 506 | | |
| 33 | | | | 519 | | |
| 34 | | | | 531 | | |
| 35 | | | | 541 | | |
| 36 | | | | 547 | | |
| 37 | | | | 550 | | |

Troubleshooting

The following table is arranged in a Trouble/Probable Cause/Remedy format. The probable causes for specific problems are listed in descending order of significance. That is, check out the first cause listed before proceeding to the next. Perform all procedures in order listed and exactly as stated. If you need additional assistance please contact Chart Technical Service at 1-800-400-4683.

| Trouble | Probable Cause | Remedy |
|--|--|--|
| Ultra-Helium Dewar container will not fill | Vent valve is closed | Open vent valve |
| with liquid | Dewar has not completely cooled down | Continue to fill |
| | Transfer line has lost vacuum (heavy frost on line) | Replace transfer line |
| | Transfer line is blocked | Purge the transfer line with helium gas |
| | Dewar has lost vacuum (heavy frost on the outside of the tank) | Return for repairs |
| | Residual LN ₂ is freezing solid (heavy tare weight noticed) | Warm container and purge with helium gas |
| Ultra-Helium Dewar container pressure | Vent valve is closed | Open the valve |
| is high and it will not vent | Relief valve has malfunctioned | Bench test the relief valve and replace if necessary |
| | Ice block has formed between the neck tubes | Remove the helium and warm up the tank |
| Ultra-Helium Dewar container has high | Shield flow valve is not working | Open the isolation valve |
| NER | | Remove and test the relief valve. Replace if necessary |
| | Dewar has been overfilled | Fill by weight to the prescribed volume |
| | Thermal-accoustical oscillation if occurring | Make sure that no optional equipment has been attached to the plumbing |
| | Dewar has marginal vacuum | Call Chart Technical Service |



Specifications

| Ultra-Helium Dewar Container Specifications | | | | | | | |
|---|-------------|--------------|--------------|--------------|---------------|--|--|
| Model | | 60 | 100 | 250 | 500 | | |
| | Part Number | 10533409 | 10533417 | 9923629 | 11202581 | | |
| Capacity | | | | | | | |
| Liquid (Gross) | liters | 66 | 110 | 275 | 550 | | |
| Liquid (Net) | liters | 60 | 100 | 250 | 500 | | |
| Performance | | | | | | | |
| NER ⁽¹⁾ | % per day | 1.75 | 1.25 | 1.0 | 1.0 | | |
| MAWP | psig / barg | 10.0 / .07 | 10.0 / 0.7 | 10.0 / 0.7 | 10.0 / 0.7 | | |
| Dimensions & Pressure R | atings | | | | ^ ^ | | |
| Relief Valve Setting | psig / barg | 10.0 / 0.7 | 10.0 / 0.7 | 10.0 / 0.7 | 10.0 / 0.7 | | |
| Secondary Relief Valve | psig / barg | 12.0 / 0.8 | 12.0 / 0.8 | 12.0 / 0.8 | 12.0 / 0.8 | | |
| Diameter | in / cm | 26 / 66 | 26 / 66 | 32.0 / 81.3 | 42.0 / 106.7 | | |
| Height | in / cm | 49.5 / 125.7 | 56.5 / 143.5 | 67.4 / 171.1 | 67.25 / 170.8 | | |
| Dip Tube Length ⁽²⁾ | in / cm | 32.5 / 82.6 | 39.5 / 100.3 | 54.4 / 138.2 | 51.5 / 98.1 | | |
| Tare Weight | lb / kg | 184 / 83 | 212 / 96 | 348 / 158 | 470 / 213.2 | | |

(1) Based on gross capacity(2) The dip tube length is measured from the tank flange to the bottom of the inner vessel

Contact Customer Service for a full line of Helium Accessories

Drawings and Parts Lists

Ultra-Helium Dewar 60 / 100





60 / 100 Parts List

| ltem | Part No. | Qty | Spares* | Description | |
|------|----------|-----|---------|---|--|
| 26 | 2013869 | 1 | 1 | Pressure Gauge - (30" Hg / 0-15 psi) | |
| 27 | 1210212 | 1 | | Coupling - 1/4" FPT | |
| 28 | 1717702 | 1 | | Ball Valve - 1/4" FPT | |
| 29 | 11366278 | 1 | 1 | Relief Valve - 1/4" MPT (1 psi) (CE if applicable) | |
| 30 | 11044869 | 1 | | Brass Tee - 1/2" FPT | |
| 31 | 11366251 | 1 | 1 | Relief Valve - 1/2" MPT (10 psi) (CE if applicable) | |
| 32 | 1312502 | 1 | | Brass Nipple - 1/2" x 3" | |
| 33 | 1717692 | 1 | | Ball Valve - 1/2" FPT | |
| 34 | 1310102 | 1 | | Brass Nipple - 1/2" x 1-1/2" | |
| 35 | 8505601 | 1 | | Top Flange Assembly | |
| 36 | 3411601 | 1 | | Vee Band Clamp | |
| 37 | 2300159 | 1 | | O-ring (Square Cut) | |
| 38 | 11491811 | 1 | 1 | Relief Valve - 1/2" MPT (12 psi) | |
| 39 | 10586201 | 1 | | Ball Valve 1/2" MPT (Fill) (PI if applicable) | |
| 40 | 2210372 | 1 | 1 | Fill Coupling Tube - 1/2" | |
| 41 | 2210072 | 1 | 1 | Fill Coupling Tube - 3/8" | |
| 42 | 2210212 | 1 | 1 | Fill Coupling Plug - 3/8" | |
| 43 | 11484282 | 2 | | Caster - 4" Dia. (Fixed) | |
| 44 | 11484303 | 2 | | Caster - 4" Dia. (Swivel) | |

*Recommended spare parts.

Ultra-Helium Dewar 250 / 500





250 / 500 Parts List

| ltem | Part No. | Qty | Spares* | Description | |
|------|----------|-------|---------|---|--|
| 51 | 2013869 | 1 | 1 | Pressure Gauge - (30" Hg / 0-15 psi) | |
| 52 | 1210212 | 1 | | Coupling - 1/4" FPT | |
| 53 | 1717702 | 1 | | Ball Valve - 1/4" FPT | |
| 54 | 11366278 | 1 | 1 | Shield Flow Valve - 1/4" MPT (1 psi) | |
| 55 | 11044869 | 1 | | Brass Tee - 1/2" FPT | |
| 56 | 11366251 | 1 | | Relief Valve - 1/2" MPT (10 psi) (CE if applicable) | |
| 57 | 1312502 | 1 | | Brass Nipple - 1/2" x 3" | |
| 58 | 1717692 | 1 | | Ball Valve - 1/2" FPT | |
| 59 | 1310102 | 1 | | Brass Nipple - 1/2" x 1-1/2" | |
| 60 | 11634252 | 1 | | Top Flange Assembly (250) 5" O.D. | |
| 60 | 11528960 | 1 | | Top Flange Assembly (500) 5" O.D. | |
| 61 | 3411601 | 1 | | Vee Band Clamp | |
| 62 | 2300159 | 1 | | O-ring (Square Cut) | |
| 63 | 11491811 | 1 | 1 | Relief Valve - 1/2" MPT (12 psi) (CE if applicable) | |
| 64 | 1717682 | 1 | | Ball Valve 3/4" Solder End (PI if applicable) | |
| 65 | 2210382 | 1 | 1 | Fill Coupling - 3/4" | |
| 66 | 2210392 | 1 | 1 | Fill Coupling Tube - 5/8" | |
| 67 | 2210372 | 1 | 1 | Fill Coupling Plug - 1/2" | |
| 68 | 1213302 | 2 | 1 | Fill Plug - 1/2" | |
| 69 | 9722169 | 2 | | P.B. Control Box | |
| 70 | 11213072 | 1 | | Street Tee - 1/4" MPT | |
| 71 | 1717819 | 1 | | Ball Valve - 1/4" MPT x 1/4" FPT (P.B.) | |
| 72 | 1013532 | 1 | | Connector - 3/16" Tube x 1/4" MPT | |
| 73 | 6910823 | 2 ft. | | Tube - 3/16" O.D. | |
| 74 | 9722009 | 1 | | P.B. Probe / Sub Assembly | |
| 75 | 2210372 | 1 | | P.B. Coupling 1/2" | |
| 76 | 2914361 | 4 | | Bolt - 1/4-20 x 1" | |
| 77 | 2914071 | 4 | | Locknut 1/4-20 | |

*Recommended spare parts.

Wiring Diagram



Schematic of Pressure Building System



Pressure Builder





Note: Layout of control box may vary slightly from drawing.

| 140.000 | Part Number | | 04. | C ===* | Description |
|---------|-------------|----------|-----|---------------|---|
| ltem | 115 Volt | 230 Volt | Qty | Spares* | Description |
| 125 | 21018407 | 21018408 | 1 | | Control Box |
| 126 | 4614989 | 4614989 | 1 | | PB Probe Plug (Female) |
| 127 | 2914202 | 2914202 | 1 | | Connector - 1/8" MPT |
| 128 | 1013372 | 1013372 | 1 | | Tube Connector - 1/8" MPT x 3/16" OD |
| 129 | 1213082 | 1213082 | 1 | | Branch Tee - 1/8" NPT |
| 130 | 21018397 | 21018397 | 1 | 1 | Pressure Switch - Set at 4 psi** |
| 131 | 21018397 | 21018397 | 1 | 1 | Pressure Switch - Set at 8 psi** |
| 132 | 14302595 | 14302595 | 1 | | Switch - 3 position |
| 133 | 21018399 | 21018398 | 2 | | Indicator - Neon Green |
| 136 | 4614329 | 14813811 | 1 | | Female Socket - 3 Wire |
| 137 | 4615069 | 14813838 | 2 | | Fuse - 1/4" 1-1/4" |
| 138 | 4614359 | N/A | 1 | | Fuse Socket - 1-1/4" x 1-1/4" for 115 Volt only |

*Recommended spare parts.

**Pressure switches must be adjusted to 4 psi or 8 psi.

