Operating Manual for MVE Liquid Nitrogen Dewars

Table of Contents	Page
English	2-7
English (SI Version)	8-13
Deutsch	14-19
Español	20-24
Italiano	25-29
Français	30-35
Português	36-41
Appendix	42
Certificate of Conformance	43

 **CHART INC.**
2200 Airport Industrial Dr., Ste 500
Ball Ground, GA 30107
U.S.A
www.chartbiomed.com

Chart Australia Pty Ltd.
ABN 21 075 909 410
Sydney Business & Technology Centre
Unit 43 / 2 Railway Parade
Lidcombe, NSW 2141, Australia

Chart BioMedical Ltd.
Unit 6 – Ashville Way
Wokingham RG41 2PL
United Kingdom

EC REP

Medical Product Services
Borngasse 20
35619 Braunsfels, Germany

Ref 11624417 Rev R 11/18

CE 0459

GENERAL DESCRIPTION

The cryopreservation vessel is a double-wall, vacuum-insulated vessel made of aluminum with a fiberglass composite neck, providing the highest efficiency possible in cryogenic temperature preservation. Use the vessel for liquid nitrogen only. Liquid oxygen is not compatible with this unit and must not be stored inside the vessel.

The MVE Liquid Nitrogen Dewar is designed with consideration for safety, durability and performance. However, mishandling of the equipment, including transport or shipping units in an orientation other than upright vertical, may damage the product. In addition, if a vessel experiences a drop, hit, or blow, it can suffer immediate or premature vacuum failure.

Upon receipt of the product, examine both the vessel and packaging for any evidence of damage during shipping. Contact the carrier within the carrier's guidelines if there are signs of shipping damage. Some MVE shipping boxes carry the Transit Tested ISTA-3A certificate stamp, shown to the right, which is helpful when making a claim against the carrier, should there be damage from shipping. Watch after the first fill for any signs of vacuum loss, such as excessive frost or sweating on the outside jacket. Some frost near the top just after filling is normal.


This high quality vacuum insulated unit is compatible with the divergent temperature extremes and broad applications of cryobiology. The Life expectancy of Liquid Nitrogen Dewar is five (5) years, Cryosystems is five (5) years.

CE Products bearing the CE marking as shown comply with the requirements of Directive 93/42/EEC 0459 concerning medical devices in the EU.

SAFETY

 **WARNING:** Liquid nitrogen is extremely cold. To avoid injury by frostbite, use extreme care whenever handling liquid nitrogen, liquid nitrogen storage or transfer vessels, or any objects which have come in contact with liquid nitrogen.



-  • Leave no area of skin exposed.
-  • Always wear proper safety attire over clothing: face shield, cryogenic gloves, and cryogenic apron.
-  • Use extreme care to prevent spilling and splashing liquid nitrogen during transfer.
-  • Always keep vessel in upright position. Do not tilt, or lay the vessel on side.
-  • Immediately remove any clothing or safety attire on which liquid nitrogen has spilled.
-  • Get immediate medical attention for any frostbite injuries due to liquid nitrogen.

 **WARNING:** The venting of nitrogen vapors may deplete oxygen in the air, possibly leading to asphyxiation or even death. Do not store or use container in areas that are small and enclosed or have poor ventilation.

 **WARNING:** Do not tightly seal liquid nitrogen container or prevent nitrogen gas from escaping. Also, excessive humidity levels or exposure to rainfall could result in freezing of the cork/cover, and possible explosion.

 **WARNING:** Never use a hollow tube to measure liquid nitrogen level. This could lead to thermal injury.

 **CAUTION:** Handle the cryopreservation vessel with care.

-  • Never overfill vessels with liquid nitrogen. Liquid nitrogen should always be below the bottom of the neck tube. Overfilling the tank may cause immediate or premature vacuum failure to occur.
-  • Never ship Liquid Nitrogen Dewar on its side or upside down. This can lead to vacuum failure.

- Remove and insert inventories carefully. Do not scratch neck tube area. Scratches can cause premature vacuum failure.
- Tampering with or removing the vacuum port will destroy vacuum and void warranty.
- Never drop, hit, or blow the unit.
- Never spill liquid nitrogen on or near vacuum port.
- Never leave the vessel in an outdoor condition.
- Keep the bottom of vessel clean and away from chemicals, fertilizers, soil, and moisture.
- Do not use MVE Liquid Nitrogen Dewar for transportation.
- All performance data published for these products is based on static conditions only. Actual performance will vary upon the nature of use. Manipulation of inventories and/or accessories along with vibration will decrease the working duration of these products.

OPERATION

 **CAUTION:** Consider the value of the inventory when storing and determine if the inventory should be split between multiple containers to reduce the risk of losing valuable and/or irreplaceable products in the unlikely event that the product fails.

 **CAUTION:** Do not use dewars for longterm storage of human biological contents. Failure to follow Chart's best operating practices, as set forth in the operating manual, can result in loss of contents.

Liquid nitrogen is extremely cold. Make sure to wear proper gear before operation. Avoid spilling liquid nitrogen over the vacuum port as this can shrink the seal and allow air to leak into the vacuum space causing premature vacuum failure. To ensure maximum performance from your MVE Liquid Nitrogen Dewar simply follow the listed steps:

1. Open container that Dewar is in, open the lid, and remove cork/cover/accessories. Lift cork/cover straight up (do not twist).
2. Fill unit to desired level. Liquid level should never pass bottom of neck tube.
 - a. If you are working with a warm vessel, it is MVE's recommendation to slowly add small amount of liquid to the bottom of unit, and allow it to sit until the liquid nitrogen stops rapidly boiling to cool the unit. Position the vacuum port facing away from the operator or other personnel.
 - b. Follow established safety practices and procedures for transferring LN2.
 - c. Fill the vessel with a funnel or transfer line when possible. Transfer using LN2 hose with phase separator or pouring container using a funnel.
 - d. If you are filling your vessel from a pressurized source, make sure it is a low pressure source (22 psi or below).
3. Replace cork/cover and allow unit to cool.
 - a. If there is excessive frost or sweating on the outside vessel after the first few hours, it would indicate either a weak or no vacuum. Examine the unit carefully.
4. Weigh unit and record.
5. Place inventories into unit, replace cork/cover and lid.

SHIPPING INSTRUCTIONS

MVE Liquid Nitrogen Dewars were primarily designed to be storage containers. Using it for transport purposes could cause permanent damage to the unit. If you need to transport your inventory at cryogenic temperature, please consider MVE Vapor Shipper series.

When transferring empty units to a location, ensure the shipping orientation is upright. Shipping units in any orientation other than upright can cause permanent damage to the tank causing loss of vacuum.

GENERAL CLEANING

Do not use any petroleum based cleaning solution.

Inside of the vessel:

Any cleaning solution that does not react with aluminum, stainless steel, or G10 composite can be used in the sanitation process of a MVE Dewar. In most cases, any household detergent or mild soap solution is suitable. Other cleaners and disinfectants that can be used safely include hydrogen peroxide, chlorine/water mixture, and denatured alcohol. The generally accepted practice of using 10% chlorine bleach with 90% water solution holds as the best method for decontamination. It is important that all surfaces being sanitized are thoroughly rinsed and that all cleaner solution residues are removed after cleaning. Allow the unit to dry completely before putting into service. It is suggested that the unit is inverted to drain and dry completely.

Outside of the vessel:

Use light dampened cloth of mild soap solution.

PREVENTIVE MAINTENANCE

The owner may follow two methods of maintenance, Continuous and/or Annual.

For Continuous Maintenance (whenever dewar is filled and/or within sight):

If the tank is warm follow steps 1 - 3 under OPERATION section above but allow the unit with full liquid to stand undisturbed for 2 hours unless dewar is already cold. Verify no major frost or condensation on outside of unit. Frost or condensation on the outside of unit would indicate either a weak or no vacuum. Inspect the cork and cover for signs of damage that may affect NER. Inspect the outside of the tank for signs of shipping damage, such as large dents, especially around the neck tube area. Be vigilant and observe the tank frequently. Remove badly damaged tanks from service.

Clean up any stain, contamination or condensation on the bottom of the unit after dumping liquid out. Keep the unit in a cool, clean and dry location for storage. These procedures can help prevent corrosion on the bottom of the vessel.

For Annual Maintenance:

1. Fill unit to approximately half full. Refer to OPERATION section step 1, 2 & 3 for details.
2. Allow unit to stand for minimal 24 hours.
3. Weigh unit and record as *First Weight*, [lb].
4. Allow filled unit to sit undisturbed for another 24 (+/-0.25) hours.
 - a. Consider the accuracy and resolution of your scale to determine if additional days are required between the first and second weights in order to obtain an accurate NER. Make sure to record number of hours between first and second weight.
5. Weigh second time and record as *Second Weight*, [lb].
6. Calculate evaporation rate by using equation below. The difference between the first weight and the second weight is the daily evaporation rate in lbs. This figure roughly signifies the normal evaporation rate, or N.E.R, [Liter/Day]

$$NER = \frac{(First\ Weight - Second\ Weight) \times 13.468}{Number\ of\ Hours}$$

If there is major frost or condensation on outside of container during this time, it would indicate either a weak or no vacuum.

Refer to the static NER specifications in Appendix A at the end of this manual. Factors such as age of unit, quantity of inventory, ambient environment, shipping condition, and use of accessories, etc. can negatively affect unit NER. If you do not find your model listed here, refer to the MVE Cryopreservation Catalog.

Note: MVE recommends the customer keeps a spare tank filled with liquid nitrogen on hand for emergency use should a tank in service become damaged or lose vacuum, thus being able to save valuable contents by transferring them into the spare tank.

MEASURING LIQUID NITROGEN LEVEL

1. Always wear proper safety attire – face shield, cryogenic gloves, and cryogenic apron.
2. Use wooden or plastic dipstick for liquid level measurement inside a unit. Using a hollow tube can cause liquid nitrogen to flash through the tube and result in thermal injury.
3. Level will be indicated by frost line, which develops when dipstick is removed.
4. You could also consider to purchase a low level alarm to monitor the liquid level.

LIQUID WITHDRAWAL

1. Always wear proper safety attire – face shield, cryogenic gloves, and cryogenic apron.
2. Liquid withdrawal from units must be accomplished by pouring or utilizing a withdrawal device. Withdrawal device pressurizes to approximately 5 psi and the pressure forces liquid up the withdrawal tube out the valve.
3. Follow the operation manual of liquid withdrawal device closely to avoid injury.

MONITORING LIQUID LEVEL

MVE supplies Liquid Level Alarms (LLA) for peace of mind when storing valuable inventory in liquid nitrogen containers. Contact Customer or Technical Service for assistance in selecting the proper LLA for your product.

Limited Warranty Statement

Chart Inc., 2200 Airport Industrial Dr., Suite 500, Ball Ground, GA 30107

General Terms:

“**New Product**” – Chart Inc. (“CHART”) warrants to the original purchaser (“Purchaser”) that each new Liquid Nitrogen Dewar (collectively, the “New Products”), shall be free from defects in materials and workmanship for a period of two (2) years from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship as listed below, except as provided in the Limited Warranty Provisions below: Spectrum and CryoSystem Series, ET11, ET24, A-200, Workhorse 2003: 3 years / XC, SC, and Lab Series, ET 23-JLM: 5 years

“**Repair/Blemished Product**” – CHART warrants to the Purchaser that all repaired equipment and factory blemished equipment (collectively, the “Repair/Blemished Products”) shall be free of defects in materials and workmanship for a period of (90) days from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship for a period of (1) year from date of shipment, except as provided below.

“**Service Part**” – CHART warrants to the Purchaser that all Service Parts and Accessory Items (each a “Service Part”) are warranted for (90) days from date of shipment to be free of defects in material and workmanship, except as provided below.

Limited Warranty Provisions:

Purchaser agrees that before this limited warranty shall become effective, Purchaser shall fully inspect each New Product, Repair/Blemished Product, or Service Part within three (3) days of delivery and before such Product is put to use. Purchaser also agrees to operate the New Product, Repair/Blemished Product, or Service Part in accordance with CHART’s operating instructions and that failure to do so shall void this limited warranty. Purchaser further agrees that any claim for breach of warranty must be made in writing within 60 days of discovery of a purported defect. CHART will not be responsible for any alleged breach of warranty, which, as a result of CHART’s inspection, CHART determines to have arisen from a cause not covered by this limited warranty. In this case, CHART will charge the purchaser a nominal fee to repair the unit.

This limited warranty does not apply to: (A) Normal routine service items; (B) Repair or replacement necessitated by misuse, abuse, accident, or repairs made by persons other than CHART or persons not authorized by CHART, (C) Use of external equipment or parts with the New Product, Repair/Blemished Product, or Service Parts other than those approved by CHART, (D) Defects caused by effects of normal wear and tear; and (E) Acts of God, or other causes not within the control of CHART.

If Purchaser believes that a New Product, Repair/Blemished Product, or Service Part does not comply with the limited warranty stated above, Purchaser should contact CHART at the address stated above, describing the problem and providing proof of the date of purchase. If directed by CHART, Purchaser shall return the New Product, Repair/Blemished Product, or Service Part freight prepaid, properly packaged in a CHART approved shipping container and properly identified by a Return Material Authorization Number issued by CHART. New Products, Repair/Blemished Product, or Service Parts returned without a Return Material Authorization Number will be refused and returned at Purchaser’s expense.

The remedies available for any breach of this limited warranty are limited to repair or replacement of the defective New Product, Repair/Blemished Product, Service Part, or refund of the purchase price, at the sole discretion of CHART. CHART warrants that replacement or repaired New Product, Repair/Blemished Product, or Service Part shall be free from defects in material and workmanship for the duration of the unexpired portion of the original warranty or ninety (90) days from the date of re-shipment to Purchaser, whichever is longer. PURCHASER'S RECOVERY FROM CHART FOR ANY CLAIM SHALL NOT EXCEED PURCHASER'S PURCHASE PRICE FOR THE NEW PRODUCT, REPAIR/BLEMISHED PRODUCT, OR SERVICE PART GIVING RISE TO SUCH CLAIM, IRRESPECTIVE OF THE NATURE OF THE CLAIM, WHETHER IN CONTRACT, TORT, WARRANTY, NEGLIGENCE, STRICT LIABILITY OR OTHERWISE. CHART SHALL NOT BE LIABLE FOR AND PURCHASER SHALL INDEMNIFY, DEFEND AND HOLD CHART HARMLESS FROM ANY CLAIMS BASED ON CHART'S COMPLIANCE WITH PURCHASER'S DESIGNS, SPECIFICATIONS OR INSTRUCTIONS, OR MODIFICATION OF ANY PRODUCTS BY PARTIES OTHER THAN CHART, OR USE IN COMBINATION WITH OTHER PRODUCTS.

PURCHASER SHALL NOT IN ANY EVENT BE ENTITLED TO, AND CHART SHALL NOT BE LIABLE FOR, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY NATURE INCLUDING, WITHOUT LIMITATION, BUSINESS INTERRUPTION COSTS, REMOVAL AND/OR REINSTALLATION COSTS, REPROCUREMENT COSTS, LOSS OF PROFIT OR REVENUE, LOSS OF DATA, PROMOTIONAL OR MANUFACTURING EXPENSES, OVERHEAD, INJURY TO REPUTATION OR LOSS OF CUSTOMERS, EVEN IF CHART HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

EXCEPT FOR THIS LIMITED WARRANTY, CHART HAS MADE NO WARRANTIES OR REPRESENTATIONS, EXPRESSED OR IMPLIED, AND HEREBY DISCLAIMS ANY OTHER WARRANTIES, INCLUDING, BUT NOT LIMITED TO, IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NO REPRESENTATION OR STATEMENT OF CHART MAY CHANGE OR ALTER THIS LIMITED WARRANTY.

Any claims for breach of this limited warranty shall be governed by Georgia law and without regard to conflict of law rules thereof and must be brought in a state or federal court in Georgia.

Some states do not allow limitations on implied warranties or on incidental or consequential damages, so the above limitations may not apply. This limited warranty gives the Purchaser specific legal rights. Purchaser may also have other rights, which vary from state to state.

This product may be covered by one or more patents, US and international. Please visit our website below for the listing of applicable patents:

Pat. patents.gtls.io

Supplement to Chart Inc. Limited Warranty Statement for Sales in Australia (AI)

General Terms:

"New Product" – Chart Inc. ("CHART") warrants to the original purchaser only (meaning the party from whom CHART accepts the purchase order), whether such purchaser is a wholesaler, distributor or end user ("Purchaser"), that each new Liquid Nitrogen Dewar (collectively, the "New Products"), shall be free from defects in materials and workmanship for a period of two (2) years from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship as listed below, except as provided in the Limited Warranty Provisions below: Spectrum and CryoSystem Series, ET11, ET24, A-200, Workhorse 2003: 3 years / XC, SC, and Lab Series, ET 23-JLM: 5 years.

"Repair/Blemished Product" – CHART warrants to the Purchaser that all repaired equipment and factory blemished equipment (collectively, the "Repair/Blemished Products") shall be free of defects in materials and workmanship for a period of (90) days from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship for a period of (1) year from date of shipment, except as provided below.

"Service Part" - CHART warrants to the Purchaser that all Service Parts and Accessory Items (each a "Service Part") are warranted for (90) days from date of shipment to be free of defects in material and workmanship, except as provided below.

Sales of Product in Australia:

Where the Purchaser is the end user of the Product, the following additional information is provided in relation to the warranty against defects (the Warranty) provided by CHART (as set out in the Limited Warranty Statement).

1. In circumstances where a New Product, Repair/Blemished Product or Service Part does not comply with the Warranty, CHART will at its expense repair or replace that defective New Product, Repair/Blemished Product or Service Part, or refund the purchase price at the sole discretion of CHART. CHART warrants that the replacement or repaired New Product, Repair/Blemished Product or Service Part shall be free from defects in materials and workmanship for the duration of the unexpired portion of the Warranty, or 90 days from the date of reshipment to the purchaser, whichever is longer.
2. The Warranty is subject to the Limited Warranty Provisions set out below and all exclusions and limitations set forth in the CHART Limited Warranty Statement.
3. To claim the Warranty, the Purchaser must contact CHART's affiliate in Australia for further instruction, including the return of the defective New Product Repair/Blemished Product or Service Part (at the Purchaser's expense), to the following address:

Chart Australia Pty Ltd.
Unit 43/2 Railway Parade
Lidcombe, NSW 2141, Australia
Attention: Customer Service Department
Phone number: (02) 9749-4333
Email address: customerservice.australia@chartindustries.com

4. The Purchaser must otherwise comply with the requirements set out in the Limited Warranty Statement when making a claim under the Warranty.
5. Notwithstanding the Limited Warranty Statement, our Products come with guarantees that cannot be excluded under the Australian Consumer Law (ACL). You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the Products repaired or replaced if the Products fail to be of acceptable quality and the failure does not amount to a major failure.
6. The provider of the Warranty is CHART. Relevant contact details are set out below for CHART (to claim the Warranty, please use the contact details provided in paragraph 3 above):
 - (a) Customer Service Department
 - (b) Business Address:
Chart Inc.
2200 Airport Industrial Drive, Suite 500
Ball Ground, GA 30107 (USA)
 - (c) Phone number: 770-721-7759
 - (d) Email address: customerservice.usa@chart-ind.com

Repair of Products Sold in Australia

Where the Purchaser is the end user of the Product, the following additional information is provided in relation to the repair of Products by CHART:

The Products may contain "user generated data" as defined under the Australian Consumer Law (ACL). Repair of Products may result in the loss of user generated data.

This product may be covered by one or more patents, US and international. Please visit our website below for the listing of applicable patents:

Pat. patents.gtls.io

Operating Manual for MVE Liquid Nitrogen Dewars (SI Version)

M.D.D. Representative: Medical Product Services, Borngasse 20, 35619 Braunsfels, Germany

GENERAL DESCRIPTION

The cryopreservation vessel is a double-wall, vacuum-insulated vessel made of aluminum with a fiberglass composite neck, providing the highest efficiency possible in cryogenic temperature preservation. Use the vessel for liquid nitrogen only. Liquid oxygen is not compatible with this unit and must not be stored inside the vessel.

The MVE Liquid Nitrogen Dewar is designed with consideration for safety, durability and performance. However, mishandling of the equipment, including transport or shipping units in an orientation other than upright vertical, may damage the product. In addition, if a vessel experiences a drop, hit, or blow, it can suffer immediate or premature vacuum failure.

Upon receipt of the product, examine both the vessel and packaging for any evidence of damage during shipping. Contact the carrier within the carrier's guidelines if there are signs of shipping damage. Some MVE shipping boxes carry the Transit Tested ISTA-3A certificate stamp, shown to the right, which is helpful when making a claim against the carrier, should there be damage from shipping. Watch after the first fill for any signs of vacuum loss, such as excessive frost or sweating on the outside jacket. Some frost near the top just after filling is normal.

This high quality vacuum insulated unit is compatible with the divergent temperature extremes and broad applications of cryobiology. The life expectancy of Liquid Nitrogen Dewar is five (5) years, Cryosystems is five (5) years.

CE Products bearing the CE marking as shown comply with the requirements of Directive 93/42/EEC 0459 concerning medical devices in the EU.

SAFETY

 **WARNING:** Liquid nitrogen is extremely cold. To avoid injury by frostbite, use extreme care whenever handling liquid nitrogen, liquid nitrogen storage or transfer vessels, or any objects which have come in contact with liquid nitrogen.

- 
- Leave no area of skin exposed.
- Always wear proper safety attire over clothing: face shield, cryogenic gloves, and cryogenic apron.
- Use extreme care to prevent spilling and splashing liquid nitrogen during transfer.
- Always keep vessel in upright position. Do not tilt, or lay the vessel on side.
- Immediately remove any clothing or safety attire on which liquid nitrogen has spilled.
- Get immediate medical attention for any frostbite injuries due to liquid nitrogen.

 **WARNING:** The venting of nitrogen vapors may deplete oxygen in the air, possibly leading to asphyxiation or even death. Do not store or use container in areas that are small and enclosed or have poor ventilation.

 **WARNING:** Do not tightly seal liquid nitrogen container or prevent nitrogen gas from escaping. Also, excessive humidity levels or exposure to rainfall could result in freezing of the cork/cover, and possible explosion.

 **WARNING:** Never use a hollow tube to measure liquid nitrogen level. This could lead to thermal injury.

 **CAUTION:** Handle the cryopreservation vessel with care.

- Never overfill vessels with liquid nitrogen. Liquid nitrogen should always be below the bottom of the neck tube. Overfilling the tank may cause immediate or premature vacuum failure to occur.

- Never ship Liquid Nitrogen Dewar on its side or upside down. This can lead to vacuum failure.
- Remove and insert inventories carefully. Do not scratch neck tube area. Scratches can cause premature vacuum failure.
- Tampering with or removing the vacuum port will destroy vacuum and void warranty.
- Never drop, hit, or blow the unit.
- Never spill liquid nitrogen on or near vacuum port.
- Never leave the vessel in an outdoor condition.
- Keep the bottom of vessel clean and away from chemicals, fertilizers, soil, and moisture.
- Do not use MVE Liquid Nitrogen Dewar for transportation.
- All performance data published for these products is based on static conditions only. Actual performance will vary upon the nature of use. Manipulation of inventories and/or accessories along with vibration will decrease the working duration of these products.

OPERATION

 **CAUTION:** Consider the value of the inventory when storing and determine if the inventory should be split between multiple containers to reduce the risk of losing valuable and/or irreplaceable products in the unlikely event that the product fails.

 **CAUTION:** Do not use dewars for longterm storage of human biological contents. Failure to follow Chart's best operating practices, as set forth in the operating manual, can result in loss of contents.

Liquid nitrogen is extremely cold. Make sure to wear proper gear before operation. Avoid spilling liquid nitrogen over the vacuum port as this can shrink the seal and allow air to leak into the vacuum space causing premature vacuum failure. To ensure maximum performance from your MVE Liquid Nitrogen Dewar simply follow the listed steps:

1. Open container that Dewar is in, open the lid, and remove cork/cover/accessories. Lift cork/cover straight up (do not twist).
2. Fill unit to desired level. Liquid level should never pass bottom of neck tube.
 - a. If you are working with a warm vessel, it is MVE's recommendation to slowly add small amount of liquid to the bottom of unit, and allow it to sit until the liquid nitrogen stops rapidly boiling to cool the unit. Position the vacuum port facing away from the operator or other personnel.
 - b. Follow established safety practices and procedures for transferring LN2.
 - c. Fill the vessel with a funnel or transfer line when possible. Transfer using LN2 hose with phase separator or pouring container using a funnel.
 - d. If you are filling your vessel from a pressurized source, make sure it is a low pressure source (1.52 bar or below).
3. Replace cork/cover and allow unit to cool.
 - a. If there is excessive frost or sweating on the outside vessel after the first few hours, it would indicate either a weak or no vacuum. Examine the unit carefully.
4. Weigh unit and record.
5. Place inventories into unit, replace cork/cover and lid.

SHIPPING INSTRUCTIONS

MVE Liquid Nitrogen Dewars were primarily designed to be storage containers. Using it for transport purposes could cause permanent damage to the unit. If you need to transport your inventory at cryogenic temperature, please consider MVE Vapor Shipper series.

When transferring empty units to a location, ensure the shipping orientation is upright. Shipping units in any orientation other than upright can cause permanent damage to the tank causing loss of vacuum.

GENERAL CLEANING

Do not use any petroleum based cleaning solution.

Inside of the vessel:

Any cleaning solution that does not react with aluminum, stainless steel, or G10 composite can be used in the sanitation process of a MVE Dewar. In most cases, any household detergent or mild soap solution is suitable. Other cleaners and disinfectants that can be used safely include hydrogen peroxide, chlorine/water mixture, and denatured alcohol. The generally accepted practice of using 10% chlorine bleach with 90% water solution holds as the best method for decontamination. It is important that all surfaces being sanitized are thoroughly rinsed and that all cleaner solution residues are removed after cleaning. Allow the unit to dry completely before putting into service. It is suggested that the unit is inverted to drain and dry completely.

Outside of the vessel:

Use light dampened cloth of mild soap solution.

PREVENTIVE MAINTENANCE

The owner may follow two methods of maintenance, Continuous and/or Annual.

For Continuous Maintenance (whenever dewar is filled and/or within sight):

If the tank is warm follow steps 1 - 3 under OPERATION section above but allow the unit with full liquid to stand undisturbed for 2 hours unless dewar is already cold. Verify no major frost or condensation on outside of unit. Frost or condensation on the outside of unit would indicate either a weak or no vacuum. Inspect the cork and cover for signs of damage that may affect NER. Inspect the outside of the tank for signs of shipping damage, such as large dents, especially around the neck tube area. Be vigilant and observe the tank frequently. Remove badly damaged tanks from service.

Clean up any stain, contamination or condensation on the bottom of the unit after dumping liquid out. Keep the unit in a cool, clean and dry location for storage. These procedures can help prevent corrosion on the bottom of the vessel.

For Annual Maintenance:

1. Fill unit to approximately half full. Refer to OPERATION section step 1, 2 & 3 for details.
2. Allow unit to stand for minimal 24 hours.
3. Weigh unit and record as *First Weight*, [kg].
4. Allow filled unit to sit undisturbed for another 24 (+/-0.25) hours.
 - a. Consider the accuracy and resolution of your scale to determine if additional days are required between the first and second weights in order to obtain an accurate NER. Make sure to record number of hours between first and second weight.
5. Weigh second time and record as *Second Weight*, [kg].
6. Calculate evaporation rate by using equation below. The difference between the first weight and the second weight is the daily evaporation rate in kg. This figure roughly signifies the normal evaporation rate, or N.E.R., [Liter/Day]

$$NER = \frac{(First\ Weight - Second\ Weight) \times 29.6919}{Number\ of\ Hours}$$

If there is major frost or condensation on outside of container during this time, it would indicate either a weak or no vacuum.

Refer to the static NER specifications in Appendix A at the end of this manual. Factors such as age of unit, quantity of inventory, ambient environment, shipping condition, and use of accessories, etc. can negatively affect unit NER. If you do not find your model listed here, refer to the MVE Cryopreservation Catalog.

Note: MVE recommends the customer keeps a spare tank filled with liquid nitrogen on hand for emergency use should a tank in service become damaged or lose vacuum, thus being able to save valuable contents by transferring them into the spare tank.

MEASURING LIQUID NITROGEN LEVEL

1. Always wear proper safety attire – face shield, cryogenic gloves, and cryogenic apron.
2. Use wooden or plastic dipstick for liquid level measurement inside a unit. Using a hollow tube can cause liquid nitrogen to flash through the tube and result in thermal injury.
3. Level will be indicated by frost line, which develops when dipstick is removed.
4. You could also consider to purchase a low level alarm to monitor the liquid level.

LIQUID WITHDRAWAL

1. Always wear proper safety attire – face shield, cryogenic gloves, and cryogenic apron.
2. Liquid withdrawal from units must be accomplished by pouring or utilizing a withdrawal device. Withdrawal device pressurizes to approximately 34.5 kPa and the pressure forces liquid up the withdrawal tube out the valve.
3. Follow the operation manual of liquid withdrawal device closely to avoid injury.

MONITORING LIQUID LEVEL

MVE supplies Liquid Level Alarms (LLA) for peace of mind when storing valuable inventory in liquid nitrogen containers. Contact Customer or Technical Service for assistance in selecting the proper LLA for your product.

Limited Warranty Statement

Chart Inc., 2200 Airport Industrial Dr., Suite 500, Ball Ground, GA 30107

General Terms:

“**New Product**” – Chart Inc. (“CHART”) warrants to the original purchaser (“Purchaser”) that each new Liquid Nitrogen Dewar (collectively, the “New Products”), shall be free from defects in materials and workmanship for a period of two (2) years from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship as listed below, except as provided in the Limited Warranty Provisions below:

Spectrum and CryoSystem Series, ET11, ET24, A-200, Workhorse 2003: 3 years / XC, SC, and Lab Series, ET 23-JLM: 5 years

“**Repair/Blemished Product**” – CHART warrants to the Purchaser that all repaired equipment and factory blemished equipment (collectively, the “Repair/Blemished Products”) shall be free of defects in materials and workmanship for a period of (90) days from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship for a period of (1) year from date of shipment, except as provided below.

“**Service Part**” – CHART warrants to the Purchaser that all Service Parts and Accessory Items (each a “Service Part”) are warranted for (90) days from date of shipment to be free of defects in material and workmanship, except as provided below.

Limited Warranty Provisions:

Purchaser agrees that before this limited warranty shall become effective, Purchaser shall fully inspect each New Product, Repair/Blemished Product, or Service Part within three (3) days of delivery and before such Product is put to use. Purchaser also agrees to operate the New Product, Repair/Blemished Product, or Service Part in accordance with CHART’s operating instructions and that failure to do so shall void this limited warranty. Purchaser further agrees that any claim for breach of warranty must be made in writing within

60 days of discovery of a purported defect. CHART will not be responsible for any alleged breach of warranty, which, as a result of CHART’s inspection, CHART determines to have arisen from a cause not covered by this limited warranty. In this case, CHART will charge the purchaser a nominal fee to repair the unit.

This limited warranty does not apply to: (A) Normal routine service items; (B) Repair or replacement necessitated by misuse, abuse, accident, or repairs made by persons other than CHART or persons not authorized by CHART, (C) Use of external equipment or parts with the New Product, Repair/Blemished Product, or Service Parts other than those approved by CHART, (D) Defects caused by effects of normal wear and tear; and (E) Acts of God, or other causes not within the control of CHART.

If Purchaser believes that a New Product, Repair/Blemished Product, or Service Part does not comply with the limited warranty stated above, Purchaser should contact CHART at the address stated above, describing the problem and providing proof of the date of purchase. If directed by CHART, Purchaser shall return the New Product, Repair/Blemished Product, or Service Part freight prepaid, properly packaged in a CHART approved shipping container and properly identified by a Return Material Authorization Number issued by CHART. New Products, Repair/Blemished Product, or Service Parts returned without a Return Material Authorization Number will be refused and returned at Purchaser’s expense.

The remedies available for any breach of this limited warranty are limited to repair or replacement of the defective New Product, Repair/Blemished Product, Service Part, or refund of the purchase price, at the sole discretion of CHART. CHART warrants that replacement or repaired New Product, Repair/Blemished Product, or Service Part shall be free from defects in material and workmanship for the duration of the unexpired portion of the original warranty or ninety (90) days from the date of re-shipment to Purchaser, whichever is longer. PURCHASER'S RECOVERY FROM CHART FOR ANY CLAIM SHALL NOT EXCEED PURCHASER'S PURCHASE PRICE FOR THE NEW PRODUCT, REPAIR/BLEMISHED PRODUCT, OR SERVICE PART GIVING RISE TO SUCH CLAIM, IRRESPECTIVE OF THE NATURE OF THE CLAIM, WHETHER IN CONTRACT, TORT, WARRANTY, NEGLIGENCE, STRICT LIABILITY OR OTHERWISE. CHART SHALL NOT BE LIABLE FOR AND PURCHASER SHALL INDEMNIFY, DEFEND AND HOLD CHART HARMLESS FROM ANY CLAIMS BASED ON CHART'S COMPLIANCE WITH PURCHASER'S DESIGNS, SPECIFICATIONS OR INSTRUCTIONS, OR MODIFICATION OF ANY PRODUCTS BY PARTIES OTHER THAN CHART, OR USE IN COMBINATION WITH OTHER PRODUCTS.

PURCHASER SHALL NOT IN ANY EVENT BE ENTITLED TO, AND CHART SHALL NOT BE LIABLE FOR, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY NATURE INCLUDING, WITHOUT LIMITATION, BUSINESS INTERRUPTION COSTS, REMOVAL AND/OR REINSTALLATION COSTS, REPROCUREMENT COSTS, LOSS OF PROFIT OR REVENUE, LOSS OF DATA, PROMOTIONAL OR MANUFACTURING EXPENSES, OVERHEAD, INJURY TO REPUTATION OR LOSS OF CUSTOMERS, EVEN IF CHART HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

EXCEPT FOR THIS LIMITED WARRANTY, CHART HAS MADE NO WARRANTIES OR REPRESENTATIONS, EXPRESSED OR IMPLIED, AND HEREBY DISCLAIMS ANY OTHER WARRANTIES, INCLUDING, BUT NOT LIMITED TO, IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NO REPRESENTATION OR STATEMENT OF CHART MAY CHANGE OR ALTER THIS LIMITED WARRANTY.

Any claims for breach of this limited warranty shall be governed by Georgia law and without regard to conflict of law rules thereof and must be brought in a state or federal court in Georgia.

Some states do not allow limitations on implied warranties or on incidental or consequential damages, so the above limitations may not apply. This limited warranty gives the Purchaser specific legal rights. Purchaser may also have other rights, which vary from state to state.

This product may be covered by one or more patents, US and international. Please visit our website below for the listing of applicable patents:

Pat. patents.gtls.io

Supplement to Chart Inc. Limited Warranty Statement for Sales in Australia (AI)

General Terms:

"New Product" – Chart Inc. ("CHART") warrants to the original purchaser only (meaning the party from whom CHART accepts the purchase order), whether such purchaser is a wholesaler, distributor or end user ("Purchaser"), that each new Liquid Nitrogen Dewar (collectively, the "New Products"), shall be free from defects in materials and workmanship for a period of two (2) years from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship as listed below, except as provided in the Limited Warranty Provisions below: Spectrum and CryoSystem Series, ET11, ET24, A-200, Workhorse 2003: 3 years / XC, SC, and Lab Series, ET 23-JLM: 5 years.

"Repair/Blemished Product" – CHART warrants to the Purchaser that all repaired equipment and factory blemished equipment (collectively, the "Repair/Blemished Products") shall be free of defects in materials and workmanship for a period of (90) days from the date of shipment, except as provided below. CHART warrants that the dewar vacuum integrity shall be free from defects in materials and workmanship for a period of (1) year from date of shipment, except as provided below.

"Service Part" - CHART warrants to the Purchaser that all Service Parts and Accessory Items (each a "Service Part") are warranted for (90) days from date of shipment to be free of defects in material and workmanship, except as provided below.

Sales of Product in Australia:

Where the Purchaser is the end user of the Product, the following additional information is provided in relation to the warranty against defects (the Warranty) provided by CHART (as set out in the Limited Warranty Statement).

1. In circumstances where a New Product, Repair/Blemished Product or Service Part does not comply with the Warranty, CHART will at its expense repair or replace that defective New Product, Repair/Blemished Product or Service Part, or refund the purchase price at the sole discretion of CHART. CHART warrants that the replacement or repaired New Product, Repair/Blemished Product or Service Part shall be free from defects in materials and workmanship for the duration of the unexpired portion of the Warranty, or 90 days from the date of reshipment to the purchaser, whichever is longer.
2. The Warranty is subject to the Limited Warranty Provisions set out below and all exclusions and limitations set forth in the CHART Limited Warranty Statement.
3. To claim the Warranty, the Purchaser must contact CHART's affiliate in Australia for further instruction, including the return of the defective New Product Repair/Blemished Product or Service Part (at the Purchaser's expense), to the following address:

Chart Australia Pty Ltd.
Unit 43/2 Railway Parade
Lidcombe, NSW 2141, Australia
Attention: Customer Service Department
Phone number: (02) 9749-4333
Email address: customerservice.australia@chartindustries.com

4. The Purchaser must otherwise comply with the requirements set out in the Limited Warranty Statement when making a claim under the Warranty.
5. Notwithstanding the Limited Warranty Statement, our Products come with guarantees that cannot be excluded under the Australian Consumer Law (ACL). You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the Products repaired or replaced if the Products fail to be of acceptable quality and the failure does not amount to a major failure.
6. The provider of the Warranty is CHART. Relevant contact details are set out below for CHART (to claim the Warranty, please use the contact details provided in paragraph 3 above):
 - (a) Customer Service Department
 - (b) Business Address:
Chart Inc.
2200 Airport Industrial Drive, Suite 500
Ball Ground, GA 30107 (USA)
 - (c) Phone number: 770-721-7759
 - (d) Email address: customerservice.usa@chart-ind.com

Repair of Products Sold in Australia

Where the Purchaser is the end user of the Product, the following additional information is provided in relation to the repair of Products by CHART:

The Products may contain "user generated data" as defined under the Australian Consumer Law (ACL). Repair of Products may result in the loss of user generated data.

This product may be covered by one or more patents, US and international. Please visit our website below for the listing of applicable patents:

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Medical Product Services, Borngasse 20, 35619 Braunsfels, Germany

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