

Hi, just a reminder that you're receiving this email because you have expressed an interest in Chart Industries. Don't forget to add taylor.erwin@chartindustries.com to your address book so we'll be sure to land in your inbox!

You may [unsubscribe](#) if you no longer wish to receive our emails.



MVE Chart Tech Tips

In This Issue

[Product Information](#)
[Helpful Hints / FAQs](#)
[Accessories](#)

PDF Version

To view Tech Tips in PDF,
[click here.](#)

Join Our List

[Join Our Mailing List!](#)

December 2013

PRODUCT INFORMATION

MVE Chart LN2 Tx

The LN2 Tx is designed specifically to remove excess liquid nitrogen from one tank to another tank safely and easily. Instead of lifting tanks filled with LN2 and pouring from one tank to another, one may use the LN2 Tx to transfer the LN2. By using a pressurized nitrogen gas source, the LN2 Tx efficiently transfers liquid nitrogen without straining your back or causing structural damage from spills.

To order, use **PN 15098872**.

Note: Operating pressure of gaseous nitrogen source must be 22-50 psi (1.5-3.4 bar).



LN2 Tx Hose

Specifications

Rigid Aluminum Tube Length ft (m)	4 (1.2)
Flex Hose Length ft (m)	8 (2.4)
Overall Weight lbs (kg)	8.5 (3.9)
Minimum Neck Opening in (mm)	3.5 (89)

Features:

- Safely, easily, quickly, and efficiently transfer LN2 from one open dewar to another
- Avoid having to lift and pour out heavy dry vapor shippers
- Can be used to empty LN2 freezers for maintenance, transport, cleaning, etc.
- Eliminates need to lift and pour dewar, reducing the chance of lower back strain
- Minimal transfer losses
- Reduce long-term LN2 usage - allows you to reuse by transferring between dewars
- Approximate flow rate up to 2 LPM


Note: One must use a source of compressed nitrogen from a cylinder or piping. Compressed air may not be used.

HELPFUL HINTS / FAQs

3-Tube Temperature Sensor Assembly

MVE Chart manufactures the 3-tube temperature sensor assembly. The sensor assembly is standard on all MVE High Efficiency and Vapor Series Freezers. Older freezers that use the dual tube assembly may use the 3-tube assembly as a replacement.

The 3-tube assembly is designed to house a third party temperature sensor. As validation and redundancy is becoming more vital to our customers, this change will make it easier to accommodate their needs. The assembly will come with a removable plug installed in the third hole.



Specifications: The injection molded XENOY material will provide thermal performance as well as lasting durability. Above are pictures of the new assembly and the inner diameters of the sensor tubes.

Service Part Numbers:

Part Number	Description	Freezer Models
14248744	3-Tube Sensor Assembly - 26" Long	MVE 1842P-150F

	(660.4 mm)	
14248816	3-Tube Sensor Assembly - 39" Long (990.6 mm)	MVE 800 and 1500 Series
14248752	3-Tube Sensor Assembly - 44" Long (1117,6 mm)	MVE 1800 Series

Disabling and Enabling the TEC 2000 Audible Alarm

Q: How does one disable the audible alarm on the TEC 2000?

A: The audible alarm can be disabled by pressing and holding the [Alarm Mute] key until the appropriate prompt is displayed. Next, use the [Down Arrow] key to toggle the arrow until it is positioned to the left of "Off" and press the [Enter] key to mute the audible alarm. To enable the audible alarm, press the [Down Arrow] key to toggle to "On" and press [Enter].

TEC 2000 Temperature Calibration

Q: How does one perform the temperature calibration on the TEC 2000?

A: There are two methods to perform the temperature calibration. The first method is the "Low Temperature Range" and the second method is the "Full Range" method. This edition of Technical Tips will cover the "Low Temperature Range" calibration. For the "Full Range" instructions, please reference page 89 in the TEC 2000 Manual. The TEC 2000 Manual may be accessed [here](#).

The "Low Temperature Range" calibration method allows for sensor accuracy near +/- 1.0°C from -150°C to 199°C. This method requires the temperature probe to be submerged in LN2. Step-by-step instructions are listed below. Chart recommends that the audible alarm be disabled prior to performing this calibration (directions above).

1. Enter the "Maintenance Menu" by pressing the following five keys in the order shown: [Liquid Usage][Up Arrow][Down Arrow][Up Arrow][Enter]
2. Verify the saturated LN2 temperature versus altitude reading is correct for your application. Reference page 34 of the TEC 2000 Manual [here](#).
3. Verify the "Use Ice Water" setting in the "Maintenance Menu" is set to "NO".
4. Exit the "Maintenance Menu" and press the [Temp A] key to display the temperature for probe "A".
5. Submerge temperature probe "A" in LN2 liquid along with three to four inches of the sensor cable. Press and hold

the [Temp A] key until the display reads "Press Enter with Temp A at LN2".


6. Press the [Enter] key to finalize the calibration of sensor "A".
7. To perform the calibration of temperature probe "B", press the [Temp B] key to display the temperature of probe "B".
8. Submerge temperature probe "B" in LN2 liquid along with three to four inches of the sensor cable.
9. Press and hold the [Temp B] key until the display reads: "Press Enter with Temp B at LN2".
10. Press the [Enter] key to finalize calibration for sensor "B".

The temperature probes are now calibrated.

ACCESSORIES

SUC Canister

MVE Chart offers rack systems that can hold single unit canisters (SUCs). An SUC can hold straws or vials on canes and can be conveniently placed into the rack system. Dimensions and part numbers for SUC canisters and racks are below:



SUC Canister: PN 9710491

SUC Canister Square Racks:

- PN 11395896: Large SUC Square Rack
 - 5.25" (133.35 mm) W x 25" (635 mm) H



- PN 11392644: Mini-SUC Square Rack
 - 3" (76.2 mm) W x 25" (635 mm) H



Upcoming Events

Trade Shows

We will be present at the following trade shows and would love to see you at our booth!

Arab Health Show

January 27-30

Dubai, UAE

Chart BioMedical-CAIRE, AirSep, and MVE products

Stand #1F38

Contact Us

For ordering information contact Customer Service

United Kingdom

+44 (0) 1189 367060

France

+33 (0) 05 61 42 94 11

Germany

+49 (0) 202 739 55429

Italy

+39 049 784 0142

customerservice.europe@chartindustries .com

Australia/New Zealand

+61 2 9749 4333

customerservice.australia@chartindustries.com

Asia/Pac Rim

770.721.7759

csasia@chartindustries.com

The Americas

770.721.7759

US Toll Free

800.482.2473

customerservice.usa@chartindustries.com

For technical information contact Technical Service:

United Kingdom	+44 (0) 1189 367060
France	+33 (0) 05 61 42 94 11
Germany	+49 (0) 202 739 55407
Italy	+39 049 644 498
	techservice.europe@chartindustries.com
Asia/Pac Rim	+61 2 9749 4333
	techservice.asia@chartindustries.com
The Americas	770.721.7759
	techservice.usa@chartindustries.com
US Toll Free	800.482.2473
	techservice.usa@chartindustries.com

[Forward email](#)



Try it FREE today.

This email was sent to snystrom@sequal.com by taylor.erwin@chartindustries.com | [Update Profile/Email Address](#) | Instant removal with [SafeUnsubscribe™](#) | [Privacy Policy](#).

Chart Industries | 2200 Airport Industrial Dr. | Ste. 500 | Ball Ground | GA | 30107