

MVE-2-35U

Product Specifications



MVE-2-35U	
Maximum Storage Capacity	
Number of Canisters	3
Number of 1/2 cc Straws 1 Level Bulk	180
Number of 1/4 cc Straws 1 Level Bulk	360
Performance	
LN2 Capacity w/o Inventory L est.	2.0
Static Evaporation Rate* L/day	0.09
Static Holding Time days	23
Unit Dimensions	
Neck Opening in. (mm)	1.4 (35)
Overall Height in. (mm)	15.55 (395)
Outer Diameter in. (mm)	7.7 (195)
Canister Height in. (mm)	4.7 (120)
Canister Diameter in. (mm)	0.98 (25)
Weight Empty lb. (kg)	7.05 (3.2)

* Static evaporation rate and static holding time are nominal. Actual rate and holding time will be affected by the nature of container use, atmospheric conditions, and manufacturing tolerances.

Conforms to MDD 93/42/EEC, the Medical Device Directive for the EU.

THREE Year Vacuum Warranty

Features include

- Engineered specifically for AI technician use
- Precision manufacturing design
- Locking cover
- Durable strap for easy carrying

The MVE-2-35U is the realization of precision robotic manufacturing, new technology and MVE's signature quality. Each component has been carefully engineered to produce a global low-cost tank while maintaining the consistent quality necessary to carry the MVE name. Statistical process controls and robotic welding and wrapping techniques allow MVE to match performance to price, while durability and function remain the same.

The MVE-2-35U is designed specifically with the elite AI technician that understands quality backed by the MVE brand and made in the USA. The lightweight profile allows for ease of transportation while ensuring quality and durability.



Visit www.chartbiomed.com for more information.



© 2020 Chart Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

Ref ML-A10046 A