



MVE 808  
BAXTER 4R9953 LAYOUT




1 DIVIDER  
80 FRAMES  
240 CASSETTES

MVE 808  
BAXTER 4R9955 LAYOUT




1 DIVIDER  
62 FRAMES  
186 CASSETTES

MVE 808  
GAMBRO DF200 LAYOUT




1 DIVIDER  
48 FRAMES  
144 CASSETTES

MVE 808  
GAMBRO DF700 LAYOUT




1 DIVIDER  
29 FRAMES  
87 CASSETTES

MVE 808  
SUC-1 LAYOUT  
WITHOUT OPTIONAL DIVIDERS



CAPACITY  
61 SUC-1 BOXES  
1,403 C-1 CANES, 8,418 1.0ml VIALS  
1,464 C-2 CANES, 8,784 1.2ml VIALS  
976 C-7 CANES, 4,880 2.0ml VIALS  
18,910 1/2cc STRAWS

MVE 808  
SUC-1 LAYOUT  
WITH OPTIONAL DIVIDERS



CAPACITY  
52 SUC-1 BOXES  
1,196 C-1 CANES, 7,176 1.0ml VIALS  
1,248 C-2 CANES, 7,488 1.2ml VIALS  
832 C-7 CANES, 4,160 2.0ml VIALS  
16,120 1/2cc STRAWS