

LCR CABINET

LIQUID CYLINDER REPAIRS CABINET

chartparts.com

Improve your liquid cylinder repair technician's performance with the LCR cabinet. Conveniently locating the right part for the right cylinder improves cylinder repair quality and safety. Integrated with the Chart Parts web site, restocking is just a click away!

Save time and repair your liquid cylinders up to OEM standards with Chart Parts and the LCR cabinet system.



LCR 1 Cabinet (shown) - 60" w x 84" h x 24" d

LCR 2 Cabinet - 36" w x 72" h x 24" d



Parts located in bins by item type.

PRODUCT HIGHLIGHTS

- Rugged industrial grade cabinet for long-life durability and security
- Cabinet has additional storage space for tools, parts and computer
- Parts organized by use or item type for fast location
- Individually labeled bins improves proper part identification necessary for safe and reliable repair
- Bin labels cross-reference to Chart Parts web site ordering form for convenient restocking
- **LCR 1:** Includes more than \$12,000 in Liquid Cylinder parts with 167 bins and a "FREE" storage cabinet as shown in the picture (\$1,899 Value)
- **LCR 2:** Includes a smaller parts cabinet with 102 bins



Innovation. Experience. Performance.®

LCR CABINET

LIQUID CYLINDER REPAIRS CABINET

chartparts.com



Toolbox storage for added convenience. (toolbox not included)

Tethered clipboard with a list of parts to track inventory.

Store laptop in cabinet for easy reordering from ChartParts.com. (laptop not included)


Small bins with labels for easy reordering.

Large bins with labels for easy reordering.

LCR 1 cabinet with doors open show organized bins with extra storage space.



Labels on each bin simplify reordering from ChartParts.com.



Lockable handle on door for security.



LCR cabinet system interfaces with Chart Parts web site for easy restocking.

Chart Inc.
U.S.: 1-800-400-4683
Worldwide: 1-952-758-4484



© 2016 Chart Inc.
P/N 13411141
www.chartparts.com