

VLCD-950 & 1900

BULK CO₂ DELIVERY SYSTEM

The VLCD bulk CO₂ delivery system is an affordable alternative for transporting liquid CO₂. Ideal for start-up installations, hot shot deliveries and remote operations, the VLCD provides easy delivery to Chart Bulk CO₂ systems. The system is mounted in a secure mobile pallet base with all the interconnecting piping and controls that are easily accessible by the driver. Heat management of the liquid CO₂ is controlled with a vacuum-insulated jacket for long hold times in periods of non-use and external pressure building systems for fast pressure recovery.

**Optional: Small accurate deliveries are performed with our exclusive FlowCom® meter system and a convenient ticket printer.*



Differential Pressure Flow Sensor & Meter



Integrated Ticket Printer



Lower Mounted Controls

PRODUCT HIGHLIGHTS

- Fills CO₂ beverage tanks efficiently from vehicle with DOT-SP 14424 exemption
- Differential Pressure Flow Sensor and Meter for accurate billing in small deliveries
- Integrated ticket printer standard for fast billing with date stamp
- Lower mounted tank controls allows easy driver access
- Palletized steel cage protects equipment and allows easy mobility with integrated fork lift slots
- Stainless steel durable tank construction with five year vacuum warranty

SYSTEM REQUIREMENTS

- 12 VDC power source
- Truck sized for full product load
- Not for use indoors unless equipped with pipeaway kit (sold separately)

Dimensions	VLCD-950		VLCD-1900	
	Length	Width	Length	Width
Length	43 in	1,092 mm	65 in	1,651 mm
Width	34 in	864 mm	46 in	1,168 mm
Height	71.25 in	1,810 mm	77 in	1,956 mm
MAWP	350 psig	24 barg	350 psig	24 barg
Capacity	950 lbs	431 kg	1,900 lbs	864 kg
Tare Weight	1,150 lbs	522 kg	2,350 lbs	1,066 kg
Full Weight	2,180 lbs	989 kg	4,250 lbs	1,932 kg