



Innovation. Experience. Performance. ®



CRYOGENIC OIL FIELD EQUIPMENT

MOBILE NITROGEN TANKS, SEMI-TRAILERS, & OFFSHORE UNITS

CRYOGENIC OIL FIELD EQUIPMENT

MOBILE NITROGEN TANKS, SEMI-TRAILERS, & OFFSHORE UNITS



Oil Field Nitrogen Tanks

Oil Field Tanks offer the best in tank design to meet customer specifications. Featuring the latest in structural designs and insulation technologies, the Oil Field Tanks are assembled with the state-of-the-art fabrication techniques to deliver maximum efficiency and durability. Available in three mounting options: saddle, truck rail and front pivot-mount to meet your specifications. Sizes range from 2,000 to 4,500 gallons and are available at any psig working pressure. Consult your Chart representative for other sizes and pressures.



Oil Field Nitrogen Semi-Trailers

Our product line of cryogenic trailers provides the maximum payloads in nitrogen service. Durability, reliability, and low maintenance are the keys to lower transportation costs. Our designs have integrated these concerns into all our models. Our engineering delivers superior light-weight construction, outstanding vacuum integrity, and thermal performance with our proprietary insulation system. Standard models range from 5,000 to 8,500 gallon capacity. Contact your Chart representative for more information on other sizes and specifications.



Offshore Oil Field Nitrogen Units

Offshore Units are used to transport bulk liquid nitrogen to offshore rigs and other industrial applications. We can engineer to your offshore unit specifications. (See photo). Constructed for rigorous demands of the oil industry, they provide a superior method of transporting and storing liquid nitrogen. These units are 20 feet long and are available in working pressures of 45 psig and 90 psig, and payload capacities from 2,000 to 3,000 gallons.

Chart Inc.
U.S.: 1-800-400-4683
Worldwide: 1-952-758-4484



© 2009 Chart Inc.
PN 14342052
www.chart-ind.com