

# AS-N (with touchscreen) PSA Oxygen Generator Parts List

(Serial Number 2192751-1 and Higher)



AS092-7 (120 V ~) and AS092-8 (220 V ~)			
	Part Number	Description	Quantity
<b>8,000 Hour/1-Year Service Parts</b>	FI049-1	Coalescing Element (MTP-95-502)	1
	FI081-1	Pre-Filter Element (FRP-95-505)	2
<b>16,000 Hour/2-Year Service Parts</b>	FI049-1	Coalescing Element (MTP-95-502)	2
	FI081-1	Pre-Filter Element (FRP-95-505)	4
	VA435-1	Equalization Valve Rebuild Kit (C131207)	1
	VA451-1	Feed Valve Rebuild Kit (C131208)	2
	VA490-1	Waste Valve Rebuild Kit (C131622)	2
	VA495-2	1/8" Drain Valve (SC8256A004V)	1
<b>Associated Items</b>	CB130-1	Oxygen Monitoring Board	
	CM178-1	5.7" HMI Color Touchscreen	
	CM194-1	Schneider PLC ( <i>Requires SN for Program</i> )	
	CM195-1	Schneider PLC Input Card	
	FI103-1	Bacteria Filter Element	
	FI127-5	Carbon Filter Element	
	FU004-1	3.0 Amp Fuse	
	FU015-1	2.0 Amp Fuse	
	FU024-1	.5 Amp Fuse	
	FU040-1	1.0 Amp Fuse	
	FU040-3	3.0 Amp Fuse	
	GA001-1	Gauge (0 – 15 psig)	
	IL057-3	LED Lamp (White)	
	MI224-1	Transducer	
	MU002-1	1/4" Muffler	
	MU030-1	1 1/2" Muffler	
	MU052-1	2" Muffler	
	MU057-1	Muffler Element	
	PS211-1	MagnaLube	
	R0651	Wilkerson R41 (Rebuild Kit Seal)	
	R1183	Wilkerson R41 (O-Ring Kit)	
	RG009-1	1 1/2" Regulator Wilkerson R41	
	RG022-1	Regulator (0 –10 psi)	
	SI001-3*	Sieve w/297.62 lbs./barrel	
	SI004-1*	Sieve w/25 lbs.	
	VA027-3	3/4" Check Valve	
	VA429-2	1 1/2" Valve Operator (Feed)	
	VA430-1	1/8" Solenoid Valve for Air Operator (Equalization)	
	VA431-1	1 1/4" Valve Operator (Equalization)	
	VA444-1	3-way Valve for Air Operator (Feed and Waste)	
VA447-1	2 1/2" Valve Operator (Waste)		

# AS-N (with touchscreen) PSA Oxygen Generator Parts List

(Serial Number 2192751-1 and Higher)



AS092-7 (120 V ~) and AS092-8 (220 V ~)			
	Part Number	Description	Quantity
<b>Associated Items</b>	VA452-1	3-way Solenoid Valve Kit (302116)	
	WH116-8	Wire Harness (Shielded)	
<b>Supplemental Equipment</b>	AN005-1	Oxygen Analyzer (Maxtec Max O <sub>2</sub> Plus)	
	AN021-1	Oxygen Analyzer (Maxtec Handi)	
	AN075-1	Oxygen Analyzer/Sensor (Maxtec Max O <sub>2</sub> Plus)	

\*The AS-N requires approximately 2,600 pounds of molecular sieve.



© 2019 AirSep Corporation. All Rights Reserved.  
AirSep reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. Any rights expressly stated herein are reserved by us, as applicable.

**Ordering Information:**  
260 Creekside Drive  
Buffalo, NY 14228-2075 U.S.A.  
Tel: (716) 691-0202 ■ Fax: (716) 691-1255  
www.airsepccd.com ■ info@airsep.com

AS-SPL-24 B