


## FLOWCOM 3000

### Flow Metering System For Cryogenic And Other Liquids

Obtain the best solution from Flow Instruments.

The Flowcom 3000 is a new generation in cryogenic metering systems from Flow Instruments. The new meter software allows programming for a perfect fit to your application – and with Flow-Pump-Smart, your pump is always protected from premature starts. Along with the most accurate and reliable meter section on the market today, the Flowcom 3000 is the ideal metering choice for atmospheric cryogenic liquids as well as natural gas, carbon dioxide, nitrous oxide, AdBlue and others.



#### Advantages:

- No moving parts
- Lifetime warranty on metering section
- High accuracy
- Long term stability
- Easy operation
- With Flow-Pump-Smart functionality:
  - Pump operation lock-out if not properly cooled down
  - Pump operation shut-down if tank runs dry
  - Hour meter with service warning (temporarily locks out pump if ignored)

#### Features:

- PLC functionality for tailor-made application software
- Integrated type approved data logger
- Optional Bluetooth interface
- Intuitive, multilingual user interface achieved by 5.7" color display and alphanumeric keypad
- Rugged zero corrosion stainless steel enclosure
- European type approval based on MID (Measurement Instruments Directive)

