

Companion[®] 1000 & 1000T

Liquid Oxygen Portables

Reliable Oxygen Delivery

The Companion[®] 1000 Portable Liquid Oxygen Unit and Companion T High-Flow Portable Liquid Oxygen Unit offer a lightweight oxygen therapy option for patients who require continuous oxygen flow. The reliable Companion 1000 unit combines comfort and convenience in a durable, simple-to-use portable liquid oxygen system. It fills quickly and easily from any Companion, Liberator[®], or HELIOS[™] Universal Reservoir. The Companion T unit offers the same design features, plus extended liters per minute settings for higher flow rates.

Features

- Small, lightweight design allows a more active lifestyle
- Easy portability
- Companion 1000 unit has 11 standard flow settings from 1/4 to 6 LPM
- Companion T unit has built-in control valve with settings from 1/2 to 15 LPM
- Both provide approximately 7.7 hours of use at 2 LPM
- Contents indicator does not require batteries
- Recessed controls and indicators minimize damage and accidental settings changes



Call 1.800.482.2473 to order!
www.CAIREmedical.com

Companion® 1000 & 1000T

Liquid Oxygen Portables



Product Specifications

	Companion 1000	Companion 1000T
LOX capacity	3.1 lb (1.41 kg)	3.1 lb (1.41 kg)
Gaseous equivalent capacity	1058 L	1058 L
Weight, empty	5.0 lb (2,27 kg)	5.7 lb (2,59 kg)
Weight, filled	8.1 lb (3,67 kg)	8.8 lb (3,99 kg)
Height	13.8 in. (351 mm)	14.5 in (368 mm)
Typical use time at 2 LPM	7.7 hours	7.7 hours
Operating pressure	22 psi (152 kPa)	22 psi (152 kPa)
Normal evaporation rate	1.5 lb/day	1.5 lb/day
Standard flow control range	Off, .25, .5, .75, 1, 1.5, 2, 2.5, 3, 4, 5, 6 LPM	Off, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15

Ordering Information

Description	Catalog Number
-------------	----------------

Companion 1000 Portable Liquid Oxygen Unit	B-775099-00
--	-------------

Companion T High-Flow Portable Liquid Oxygen Unit	B-775532-00
---	-------------

Accessories

LOX backpack	100181138
--------------	-----------

Cart, Companion 1000	B-775295-00
----------------------	-------------

Pressure gauge, tube adapter	B-775270-00
------------------------------	-------------

Wheelchair basket, transport	B-775538-00
------------------------------	-------------

CAIRE, Inc.
2200 Airport Industrial Drive, Suite 500
Ball Ground, GA 30107
Ph 770-721-7759 • Toll-Free 1-800-482-2473
Fax 888-932-2473

Chart BioMedical, Ltd.
Unit 6, Ashville Way
Wokingham RG41 2PL, United Kingdom
Ph +44(0) 1189 367060 • Fax +44(0) 1189 799245