

## Euro-Cyl Premium Transportable Liquid Cylinders



Better convenience, quality, safety and savings for your gas supply.

The units are ideal for liquid nitrogen, oxygen, argon, CO<sub>2</sub>, N<sub>2</sub>O or LNG service.

Standard product feature includes:

- Capacity from 120 up to 1000 liters
- Pressures from 4 up to 37 bar
- Variable base type foot ring, frame or caster-case
- Pressure control device
- All products utilise Chart Vacuum Technology
- They are built to last

*Euro-Cyl liquid cylinders are the most user-friendly and reliable portable tanks on the market today.*



Transportable & Stationary storage systems are reliable replacement for high-pressure cylinders

## LIQUID CYLINDRES FOR AIR GASES





## MicroBulk Distribution System Fill-on-Site Technology



The technology includes tanker „ORCA“ and liquid cylinders „Perma-Cyl“.

**ORCA** is a fast filling cryogenic tanker with ADR and EN 13 530 compliance. Submerged pump eliminates cool down losses and reduces delivery time.

**Perma-Cyl and Perma-Cyl ZX** are superior stationary liquid cylinders

Standard product features include:

- Capacity 230 up to 3000 liters for MicroBulk
- Pressures up to 37 bar
- Novel top fill float system for single hose filling without losing product
- Super insulated permanent vacuum system provides days of holding time without losing product through boil off
- The units are ideal for liquid nitrogen, oxygen, argon, CO<sub>2</sub>, N<sub>2</sub>O or LNG service

High quality workmanship and finish with excellent reliability.



Flow meter  
Dynaflow 3000  
with MID approval



*MicroBulk System is the most economical delivery system on the market.*

*MicroBulk to your householders and workshops.*

