



The VRV Group

SCOH

Mobile Gas filling station

1000 to 12000L

Mars 2017



Products and Strategy

modular system 1000L to 12000L PED/TPED/ADR with

- **High pressure pump,**
- **control panel,**
- **Vaporizer,**
- **filling manifold,**

**skid , or protection frame or trailer mounted for 200/300 bar
for LIN/LOX/LCO2**

Mobil gas filling station



600L & 1200L (6 units in 1980)

30 years of
History



12 000L (10 units in 2013)



6 000L (6 units in 1980)



3 000L (4 units in 2016)

Products and Strategy



**3000L trailer
mounted with
pump, vaporizer
and control panel**

Products and Strategy



**12 000L skid
mounted with
pump, vaporizer
and control panel**

ISO tanks and road transport units

7000 or 12000 Liter in 20ft frame

