





Valve Disc Material Formerly CPI 140



CPI 543 is the most widely selected valve disc material for use in a broad range of operating conditions.

CPI, part of the Howden group, should be consulted for the proper design and application of its specialized products and materials. For further advice and technical support please contact CPI directly.

Typical properties	Metric	Imperial
Tensile strength at 20°C	185 MPa	26,000 psi
Elongation at 20°C (%)	3-4	3-4
Coefficient of thermal expansion	32 x 10 ⁻⁶ /°C	1.8 x10 ⁻⁵ /°F
Specific gravity	1.37	1.37
Water absorption (%)	0.9	0.9
Flexural strength	260 MPa	37,000 psi
Felxural modulus	9 x 10 ³ MPa	13 x 10 ⁵ psi
Suggested gas discharge temp. limit	150°C	300°F



Note: The values above are for reference only and are not intended for specification or quality control purposes.